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National High-tech Enterprise
Jiangxi Solid Preparation Manufacturing Equipment Industrial Design Center
Pharmaceutical Solid Preparation Equipment And Information System Supplier

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YICHUN WANSHEN PHARMACEUTICAL MACHINERY CO., LTD.

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COMPANY ADDRESS

江西宜春市经济技术开发区春风路28号
NO.28 Chunfeng Road, Yichun Economic & Technological
Development Zone, Jiangxi Province, China
www.wonsen.com.cn

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DOMESTIC BUSINESS DEPARTMENT

Tel: 0086-795-2197988
Mobile: 0086-13617058788
E-mail: xiaoshou@wonsen.com.cn

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INTERNATIONAL BUSINESS DEPARTMENT

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WONSEN ENTERPRISE

Pharmaceutical solid preparation equipment solutions

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Yichun Wanshen Pharmaceutical Machinery Co., Ltd. was founded in September 2010. Wonsen are the state-level high-tech enterprise, specializing in the research, development, production, sales and service of pharmaceutical equipment and intelligent information system. It provides the global customers with turn-key solution of solid dosage equipment and information management system to make powder, granule, capsule, tablet, etc. Wonsen takes the initiative to combine the solid preparation equipment with the Internet technology. And now Wonsen has become an important industrial base for solid preparation equipment and intelligent information system in China.

Our company covers an area of 93,000 square meters, including construction area of 60,000 square meters, with industry-leading industrial design centers and digital test experience centers built according to EU GMP standard. It has successively cooperated with Zhejiang University, Hunan University, Jiangxi University of Traditional Chinese Medicine and Nanchang University. Wonsen constantly explores the key technologies and breaks through the equipment bottlenecks in the future development of pharmaceutical solid preparations. It aims to lead domestic solid preparation equipment and Internet technology applications to catch up with the world's advanced level.

Our company takes intelligent information management and lean production site management, carries forward the enterprise spirit of "being humanistic, unrelenting, pragmatic, innovative, and keeping exceeding", and advocates the business philosophy of "pursuing high quality, creating value for customers and creating benefits for our company". Adhering to the enterprise style of diligence, pragmatism, trustworthiness and striving for success, Wonsen is committed to build a century-old enterprise group up to international leading level.

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Experienced
Professional
All-round service
Problem solver

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Wonsen offers employees a platform for creating wealth and happy life. Wonsen respects every employee and holds the concept that talents are the first wealth of our company. Oriented by performance, Wonsen helps employees realize self-value through fair and just competition. Wonsen has created a harmonious and friendly working atmosphere with humanistic enterprise management and sound corporate culture.

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
Sincerity constructs harmony
and wins honors.

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Innovation is the engine for survival and development of enterprise and is the base of corporate cultural spirit of Wonsen. Wonsen is always striving to build an innovative environment, introduce a competition system, enhance innovative conception and improve innovative ability of employees. Our enterprise pays much attention to information orientation, and applies new knowledge and advanced technology into updating product. The research and development ability has been increasingly enhanced, and Wonsen has been more vigorous with new products and technology emerging. Furthermore, innovation does not merely rest on the product level, because Wonsen knows deeply that the product is only one part, not all of customer requirements! Therefore, we have been devoted to the management system reform and service mode renovation to standardize the management system and focus on every chain and detail of projects.

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The engineers are providing more solutions while creating core technologies for the industry.

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Wonsen sales center has a contingent of professional service personnel who are well-trained teamwork players. Furthermore, we have also established organizations running in a highly efficient way according to the modern enterprise system. Facing the 21st century, our sales center will fully introduce the CS (customer satisfaction) system and sincerely provide the most excellent service for merchants and users at home and abroad with perfect management and service systems by regarding customer satisfaction as the standard of evaluating our work.

万申企业 WONSEN ENTERPRISE

销售团队协作
缔造完美品质
Good teamwork leads
to exquisite product

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Management is the soul of Wonsen,
and talents are the foundation of it.

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The perfect system, highly efficient and strict operation system are the solid foundation for sound management decision and stable business development of our enterprise. Over the years, the behavior standards and rules of employees are continuously improving, and the enterprise management is moving towards a scientific, standard and professional level. Now, our company has a contingent of well-trained employees who can meet the needs for modern enterprise management and production, laying a firm foundation for future development of our company.

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The strict workshop site management, 6S location management and various operation rules for safe production have become driving force for the development of Wonsen. Wonsen are continuously intensifying the research, development and innovation of core technology, pursuing prominence and keeping improving, ensuring the excellent performance of every single product by advanced production equipment, strict quality control and precise production and inspection procedures.

精益管理 Lean management

万申企业 WONSEN ENTERPRISE

雄厚的技术
精良的装备
Strong technology and
exquisite equipment

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
Market competition is closely associated with production strength, state-of-the-art hardware facilities as well as advanced management philosophy. Deeply understanding the importance of hardware facilities, Wonsen took the lead among its peers in introducing advanced production equipment at home and abroad and ISO9001 quality management system, so as to produce high quality products conforming to FDA, CGMP and GMP standards and fully meet the requirements for industrial modernization.

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Wonsen has accumulated a large number of talents and abundant experience after many years of development. With professional technology and exquisite quality, Wonsen has attracted many domestic and overseas customers. During the cooperation and communication with them, Wonsen has become more aware of the importance of excellent quality and considerate service. The process of cooperation is also a process of accumulating our experience of growing and learning. Wonsen are striving to be better.

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Wonsen is your most
trustworthy partner

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Faithful customers are the base for Wonsen to survive. "Being sincere, faithful, pioneering and enterprising" is the marketing concept Wonsen always abides by, and it is our goal to seek for mutual progress with customers. Over the years, by consistently adhering to diversified marketing strategies, Wonsen is actively developing marketing activities revolving around customers and executing the regional sales strategy.

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Wonsen provides top-ranking professional project installation experience and management experience. By holding the work attitude of keeping improving, Wonsen provides users with good after-sales service and meets customer requirements. Professional personnel are the guarantee for good service. Our professional personnel have practical experience in equipment use and maintenance and receive systemic trainings, mastering professional knowledge and site treatment modes.

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
Scientific, professional and elaborate
manufacturing is perfectly expressed
in every procedure.

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Application

The Closed granulation line is mainly used in the powder granules of the pharmaceutical, food, chemical and other industries, which can inhibit the dust leakage and pollution generated during the process of material transfer, posing a threat to people and the environment. The closed conveying method ensures that the conveying efficiency of granulation, drying and granulation process is more scientific and reasonable, reduces the space occupied by the equipment, and makes the production process more applicable to the GMP requirements.

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Features

The Closed granulation line fully considers the compatibility of the equipment between the systems. According to the actual production capacity and dosage form design of the user, the production capacity of the whole processes can be optimally configured to ensure the continuity and stability of production. By the industrial tablet computer, automation control system can store, collect and print the process data separately, centrally control the parameter setting and monitoring of multiple devices, reflect the working status, data statistics and fault self-diagnosis in real time, making the process more stable, fully automated and intelligent.

The equipment can be operated independently or controlled by system to realize the optimal parameterized process control of the equipment, reducing the pause and dwell time, and the intelligent cooperation of the equipment makes the granulation process optimal, thus obtaining efficient granulation yield and consistency. Automatic telescopic cleaning nozzles can be set according to user requirements, so that CIP cleaning can be reliably guaranteed.

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Work principle

- ▲Process control
Based on advanced industrial tablet computer automation control, it provides complete process control and various parameter control, including graphical operation interface, online monitoring, trend chart, recipe processing, batch management, barcode processing, CIP control, and self-diagnosis processing. The management system is fully compatible according to industry standard communication network protocol control strategy and factory-wide monitoring system.
- ▲Process analysis
FDA's PAT (Process Analysis Technology) combines the design capabilities of its equipment with process engineering expertise to integrate an online (PAT) analyzer into the system to provide a true understanding of the process and help customers achieve key product quality objective so that you can actively control the key quality characteristics of your product.
- ▲Closed control
The unique Closed control consists of integrating stable sealing technology into a complete solution for pharmaceutical solids facilities, closely monitoring the quality and GMP standards applied to the customer's site, providing customized and efficient solutions for a wide range of specifications and complex projects.
- ▲Safety and environment
In order to fully comply with national and local regulations, for emission control, as well as the recycling the solvent, the exit filter and the fully enclosed equipment can be provided as needed according to the explosion-proof and anti-pressure standards.

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Application

The machine is a process machine widely applied for solid preparation production in the pharmaceutical industry. It has such functions as mixing, granulating, etc. It has been widely used in such industries as medicine, food, chemical industry, etc.

Work principle


The materials are fed into the mixing vessel through the vacuum feeder or lifting feeder. Through rotation and push by the mixing paddle at the bottom, they first move under a fluidizing status and get sufficient mixing. Then the adhesive is injected through the pressure spray gun to turn the dry powders into wet and soft materials. Meanwhile, they are made into even wet granules under the double actions of the mixing paddle and the high-speed cutter on the side wall. After granulation, the materials are discharged to the next procedure (milling or drying).

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Features

- ▲The whole machine is made of high quality stainless steel.
- ▲It is provided with special mixing barrel and conical column of a compound structure. The corners on the inner wall are subjected to smooth transition and mirror polishing with easy cleaning and without material sticking.
- ▲The mixing and granulating procedure can be completed in one step with the materials in the closed mixing barrel.
- ▲It adopts PLC control. The mixing and granulating motors are subjected to variable frequency speed adjustment. All operations can be carried out automatically or manually.
- ▲It is provided with interlocking protection device to ensure safer operation.
- ▲It is equipped with printer to automatically print all process parameters, with the original records conforming to the actual ones.
- ▲It can be equipped with wet type milling machine and swing granulator and connect with fluid bed dryer to form a fully closed granulating production line.

技术参数 Technical parameter

项目 Item	型号 Model	SHLG-100	SHLG-150	SHLG-200	SHLG-300	SHLG-400	SHLG-500	SHLG-600	SHLG-800	SHLG-1000
生产能力 Production capacity (kg/批batch)		20-40	30-60	40-80	60-120	80-160	100-200	120-240	160-320	200-400
搅拌电机功率 Mixing motor power(KW)		11	15	22	30	30	37	45	55	55
搅拌桨转速 Mixing paddle rotation speed(rpm)		0-200	0-200	0-200	0-200	0-200	0-170	0-170	0-150	0-120
制粒电机功率 Granulating motor power(KW)		4	5.5	5.5	7.5	7.5	7.5	7.5	11	11
制粒刀转速 Granulating knife rotation speed(rpm)		0-2900	0-2900	0-2900	0-2900	0-2900	0-2900	0-2900	0-2900	0-2900
整粒功率 Miller Power(Kw)		2.2	2.2	3	4	4	5.5	5.5	5.5	5.5
整粒转速 Impeller Speed(rpm)		0-900	0-900	0-900	0-900	0-900	0-900	0-900	0-900	0-900
压缩空气耗量 Compressed air consumption(m³/min)		0.4	0.4	0.6	0.6	0.6	1	1	1	1
压缩空气压力 Compressed air pressure(MPa)		0.3-0.6	0.3-0.6	0.3-0.6	0.3-0.6	0.3-0.6	0.3-0.6	0.3-0.6	0.3-0.6	0.3-0.6
设备重量 Equipment weight(kg)		1200	1500	1800	2000	2200	2500	3000	3500	4000
外形尺寸 Outline dimensions	L	2200	2365	2580	2580	2605	2780	2780	2980	3180
	L1	1800	1800	2000	2000	2000	2200	2280	2480	2680
	W	2300	2400	2400	2600	2600	2800	2800	2900	3100
	W1	1500	1600	1600	1800	1800	2000	2000	2100	2300
	H	2600	2600	3050	3050	3050	3050	3050	3250	3350
	H1	1400	1400	1400	1400	1400	1400	1400	1600	1600
	H2	1360	1360	1360	1360	1360	1370	1370	1570	1570

注：本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.

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Application


The machine is a process machine widely applied for solid preparation production in the pharmaceutical industry. It has such functions as mixing, granulating, etc. It is widely used in such industries as medicine, food, chemical industry, etc.

Work principle

The materials are fed into the mixing vessel through the vacuum feeder or lifting feeder. Through rotation and push by the mixing paddle at the bottom, they first move under a fluidizing status and get sufficient mixing. Then the adhesive is injected through the pressure spray gun to turn the dry powders into wet and soft materials. Meanwhile, they are made into even wet granules under the double actions of the mixing paddle and the high-speed cutter on the side wall. After granulation, the materials are discharged to the next procedure (milling or drying).

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
Features

- ▲The whole machine is made of high quality stainless steel.
- ▲It is provided with special mixing barrel and conical column of a compound structure. The corners on the inner wall are subjected to smooth transition and mirror polishing with easy cleaning and without material sticking.
- ▲The mixing and granulating procedure can be completed in one step with the materials in the closed mixing barrel.
- ▲It adopts PLC control. The mixing and granulating motors are subjected to variable frequency speed adjustment. All operations can be carried out automatically or manually.
- ▲It is provided with interlocking protection device to ensure safer operation.
- ▲It is equipped with printer to automatically print all process parameters, with the original records conforming to the actual ones.

技术参数 Technical parameter

项目 Item	型号 Model	SHL-100	SHL-150	SHL-200	SHL-300	SHL-400	SHL-500	SHL-600	SHL-800
生产能力 Production capacity (kg/批batch)		20-40	30-60	40-80	60-120	80-160	100-200	120-240	160-320
搅拌电机功率 Mixing motor power(KW)		11	15	22	30	30	37	45	55
搅拌桨转速 Mixing paddle rotation speed(rpm)		0-200	0-200	0-200	0-200	0-200	0-170	0-170	0-150
制粒电机功率 Granulating motor power(KW)		4	5.5	5.5	7.5	7.5	7.5	7.5	11
制粒刀转速 Granulating knife rotation speed(rpm)		0-2900	0-2900	0-2900	0-2900	0-2900	0-2900	0-2900	0-2900
压缩空气耗量 Compressed air consumption(m³/min)		0.4	0.4	0.6	0.6	0.6	1	1	1
压缩空气压力 Compressed air pressure(MPa)		0.3-0.6	0.3-0.6	0.3-0.6	0.3-0.6	0.3-0.6	0.3-0.6	0.3-0.6	0.3-0.6
设备重量 Equipment weight(kg)		1000	1300	1600	1800	2000	2200	2700	3200
外形尺寸 Outline dimensions	L	2300	2465	2580	2680	2705	2880	2880	3080
	L1	1900	1900	2000	2100	2100	2300	2300	2580
	W	1600	1600	1800	1800	2000	2150	2150	2200
	W1	900	900	1000	1000	1200	1350	1350	1400
	H	2500	2550	2930	2930	2930	3050	3050	3250
	H1	1050	1050	1280	1280	1280	1400	1400	1600
	H2	1010	1010	1240	1240	1240	1370	1370	1570

注：本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.

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Application

The machine is mainly suitable for small batch granulation process in the laboratory. The powder material and the adhesive are mixed at high speed in the bottom mixing paddle of the wet cylinder, and then cut into uniform wet granules by the side high speed cutter. The purpose of the granules is one of the process equipment commonly used in the production of solid preparations in the pharmaceutical industry, and is widely used in the pharmaceutical, food, chemical and other industries.

Work principle


The powder materials are fed into the mixing vessel manually. Through rotation and push by the mixing paddle at the bottom, they first move under a fluidizing status and get sufficient mixing. Then the adhesive is injected through the pressure spray gun to turn the dry powders into wet and soft materials. Meanwhile, they are made into even wet granules under the double actions of the mixing paddle and the high-speed cutter on the side wall. After granulation, the materials are discharged to the next procedure (milling or drying).

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Features

- ▲The whole machine is made of high quality stainless steel and is convenient to move with castors.
- ▲The mixing and granulating procedure can be completed in one step with the materials mixed and granulated in the closed mixing barrel.
- ▲It adopts PLC control. The mixing and granulating motors are subjected to variable frequency speed adjustment. All operations can be carried out automatically or manually.
- ▲It is provided with interlocking protection device to ensure safe operation.
- ▲It is equipped with printer to automatically print all process parameters, with the original records conforming to the actual ones.

技术参数 Technical parameter

项目 Item	型号 Model	SHLS-3	SHLS-6	SHLS-10	SHLS-15	SHLS-25	SHLS-50
生产能力 Production capacity (kg/batch)(kg/批)		0.6-1.2	1.2-2.4	2-4	3-6	5-10	10-20
搅拌电机功率 Mixing motor power(KW)		1.5	1.5	1.5	4	4	4
搅拌桨转速 Mixing paddle rotation speed(rpm)		0-400	0-400	0-400	0-220	0-220	0-220
制粒电机功率 Granulating motor power(KW)		1.1	1.1	1.1	2.2	2.2	2.2
制粒刀转速 Granulating knife rotation speed(rpm)		0-2900	0-2900	0-2900	0-2900	0-2900	0-2900
压缩空气耗量 Compressed air consumption(m³/min)		0.1	0.1	0.1	0.3	0.3	0.3
压缩空气压力 Compressed air pressure(MPa)		0.3-0.6	0.3-0.6	0.3-0.6	0.3-0.6	0.3-0.6	0.3-0.6
设备重量 Equipment weight(kg)		250	270	300	400	430	450
外形尺寸 Outline dimensions	L	1400	1425	1450	1525	1550	1600
	L1	1000	1000	1000	1100	1100	1100
	W	700	700	700	750	750	750
	H	1110	1120	1140	1250	1300	1350
	H1	750	750	750	750	750	750
	H2	730	730	730	720	720	720

注：本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.

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Application

The machine is a process machine widely applied for solid preparation production in the pharmaceutical industry. It has such functions as mixing, drying, granulation functions etc. It is also widely used in such industries as medicine, chemical industry, food, etc.

Work principle

After purified and heated, the air sucked from the bottom of the fluid bed granulator by the induced draft fan, and then passes through the air flow distribution plate of the material silo, so that the material powder particles or wet granules under boiling and fluidizing status in the raw material container. After the granules are preheated by the purified and heated air and mixed, the adhesive solution is atomized and sprayed during granulating to make several particles agglomerate into granules containing the adhesive. Because the material are continuously dried by hot air, the moisture in the granules is evaporated, and the adhesive is solidified. This process will be carried out continuously and repetitively to form ideal and even micro-porous ball granules. The Fluid bed granulating method is a comprehensive application of spraying technology and fluidizing technology. It completes the traditional mixing, granulating and drying processes in a machine to realize one-step granulation.

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
Features

- ▲This machine is provided with the closed sampling function.
- ▲The air intake heating system can adopt two heating modes such as electric heating or steam heating.
- ▲One machine can be equipped with many silos.
- ▲Adopt PLC automatic control system.
- ▲Equipped with the advanced CIP cleaning system.
- ▲Meet GMP requirements.

技术参数 Technical parameter

项目 Item	型号 Model	FL-15	FL-30	FL-60	FL-120	FL-200	FL-300	FL-500
容器容积 Vessel volume(L)		45	100	220	330	577	980	1530
生产能力 Production capacity(kg/批batch)		5-15	15-30	30-60	60-120	120-200	200-300	300-500
风机功率 Fan power(KW)		7.5	11	18.5/22	22/30	30/37	37/45	75
电加热功率 Electric heating power(KW)		30	30	30	45	80	90	120
蒸汽压力 Steam pressure(MPa)		0.4-0.6	0.4-0.6	0.4-0.6	0.4-0.6	0.4-0.6	0.4-0.6	0.4-0.6
蒸汽耗量 Steam consumption(kg/h)		180	180	300	360	420	480	677
压缩空气压力 Compressed air pressure(MPa)		0.4-0.6	0.4-0.6	0.4-0.6	0.4-0.6	0.4-0.6	0.4-0.6	0.4-0.6
压缩空气耗量 Compressed air consumption(m³/min)		0.4	0.9	0.9	1.0	1.0	1.5	1.8
主机重量 Weight of main machine(kg)		800	1000	1200	1400	2000	2500	3500
外形尺寸 Outline dimensions	H	3000	3234	3953	4508	4640	5165	5812
	H1	1850	1850	1850	1850	1850	1850	1850
	ΦD	806	806	1106	1306	1306	1608	2008
	W	984	984	1340	1540	1540	1840	2240

注：本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.

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Application

The machine is a process machine widely applied for solid preparation production in the pharmaceutical industry. It has such functions as mixing, drying, etc. It is also widely used in such industries as medicine, chemical industry, food, etc.

Work principle

After purified and heated, the air sucked from the bottom of the fluid bed granulator by the induced draft fan, and then passes through the air flow distribution plate of the material silo, so that the material powder particles or wet granules under boiling and fluidizing status in the raw material container. After the granules are preheated by the purified and heated air and mixed, their moisture is evaporated. This process will be carried out continuously and repetitively to form ideal and uniform dry granules.

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

Application

This machine is one of the process equipment commonly used in the production of solid preparations in the pharmaceutical industry. It has many functions such as mixing, drying, granulating and coating. it is also widely used in the pharmaceutical, chemical, food and other industries.

Work principle

The machine is composed of a top spray system, a bottom spray system, an air treatment system, a heating system, a control system, an exhaust dust removal system and the like. The top spray is a granulation drying function, and the bottom spray is a coating function. The above functions are all realized on the same boiling granulator, and the user can select the functional requirements according to the actual production process.

Top spray principle:

After purified and heated, the air is distributed through the lower part of the fluid bed dryer through the air blower, so that the material powder particles or wet granules under boiling and fluidizing status in the raw material container. After the granules are preheated by the purified and heated air and mixed, their moisture is evaporated. This process will be carried out continuously and repetitively to form ideal and even dry granules.

Bottom spray coating principle:


When the fluidized bed is coated, the purified air process through the air distribution plate of the material silo, causing the material to be fluidized in the raw material container while atomizing the coating liquid into the package, wrapped on the surface of the material. Dry to the required process requirements.

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Features

- ▲That one machine is able to process granulation by top spray and coating by bottom spray.
- ▲This machine is provided with the closed sampling function.
- ▲The air intake heating system can adopt two heating modes such as electric heating or steam heating.
- ▲Reliable temperature control system in inlet air handling system, deviation is with ± 1.5℃.
- ▲Equipped with the advanced CIP cleaning system.
- ▲It adopted PLC automatic control system.
- ▲Meet GMP requirements .

技术参数 Technical parameter

项目 Item		型号 Model	DPL-30	DPL-60	DPL-120	DPL-200	DPL-300
容器容积 Vessel volume(L)	底喷 Bottom spray		80	150	330	417	765
	顶喷 Top spray		100	220	330	577	980
生产能力 Production capacity (kg/批batch)	底喷 Bottom spray		15-20	15-30	30-60	60-90	100-200
	顶喷 Top spray		10-30	30-60	60-120	120-200	200-300
风机功率 Fan power(KW)			11	18.5/22	22/30	30/37	37/45
电加热功率 Electric heating power(KW)			30	30	45	80	90
蒸汽压力 Steam pressure(MPa)			0.4	0.4	0.4	0.4	0.4
蒸汽耗量 Steam consumption(kg/h)			180	300	360	420	481
压缩空气压力 Compressed air pressure(MPa)			0.4-0.6	0.4-0.6	0.4-0.6	0.4-0.6	0.4-0.6
压缩空气耗量 Compressed air consumption(m³/min)			0.6	0.9	1	1	1.5
主机重量 Weight of main machine(kg)			1200	1600	1800	2300	2500
外形尺寸 Outline dimensions (mm)	H		3374	4353	4908	5040	5865
	ΦD		806	1106	1306	1306	1608
	W		984	1340	1540	1540	1840

注：本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.

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
Application

The machine is a multi-functional granulator integrating top spray, bottom spray and side spray. It is mainly suitable for small batch processing in the laboratory. According to the process requirements can be realized material mixing, drying, granulation, coating, pelleting, etc. Widely used in pharmaceutical, food, chemical and other industries.

Work principle

The machine is composed of air inlet system, fluid bed system, top spray system, bottom spray system, side spray system and exhaust and dust filter system. The top spray provides the granulating function, the bottom spray provides the pellet coating function and the side spray provides the pelletizing function. The above functions are realized in one machine. Users can select the optional functions according to actual process requirements.

技术参数 Technical parameter

项目 Item	型号 Model	DPLS-0.2	DPLS-0.5	DPLS-1	DPLS-2	DPLS-3	DPLS-5
容积 Volume(L)	顶喷 Top spray	1.25	3	5.6	5.8	9	38
	底喷 Bottom spray	1.25	3	14	16	9	22
	侧喷 Side spray	-	-	5.9	7	8.5	15
生产能力 Production capacity(kg/h)	顶喷 Top spray	0.2	0.5	1.5	2	3	3-5
	底喷 Bottom spray	0.2	0.3	1.2	2	2	1-3
	侧喷 Side spray	-	-	1	1	2	2.4
压缩空气 Compressed air	压力 Pressure(Mpa)	0.4-0.6					
	耗气量(m³/min) Air consumption	0.15	0.15	0.15	0.15	0.8	1
电加热功率 Electric heating power(KW)		2	3	4	4.5	9	12
热风调温范围 Hot air temperature adjusting range(℃)		室温~120 Room temperature to 120					
热空气过滤精度 Hot air filtering precision		亚高效或高效 Sub-efficient or effective					
风机功率 Fan power(kw)		0.75	1.75	2.2	4	4	5.5
外形尺寸 Outline dimensions	W	700	855	855	855	885	1100
	H	900	1815	2220	2275	2520	2825
	L	1000	1200	1600	1800	1800	1900
	ΦD	159	214	256	300	300	400
设备重量 Equipment weight(kg)		600	650	700	750	800	850

注：本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.

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Application

The JTZ series lifting dry mill can directly make the materials in the fluid bed granulating silo into desired granules and send them into the mixing barrel.

Work principle



The machine is composed of the fixed lifting transferring machine and milling machine. During operation, the conical hopper is assembled onto the silo. The silo is lifted, turned and the column is rotated horizontally to the required position. The silo is lowered to the work height, the milling machine is started, and with the opening of the butterfly valve, the materials are transferred to the next procedure.

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Features

- ▲The machine adopts electrically and hydraulically integrated control technology to realize stable and reliable ,safe performance. It uses the rocker operation system, which makes the use convenient and flexible.
- ▲The granules milling machine can be lowered to any work height together with the conical copper for the convenience of connection with the next procedure.
- ▲Procedure transferring is carried out under the closed status, avoiding dust and causing no pollution.
- ▲Different specification of screen can be replaced for the milling machine.
- ▲Meet GMP requirements.

技术参数 Technical parameter

项目 Item	型号 Model	JTZ-120	JTZ-200	JTZ-300	JTZ-500
净负载 Net load(kg)		120	200	300	500
提升电机功率 Lifting motor power(kw)		2.2	3.0	3.0	3.0
翻转电机功率 Turning motor power(kw)		0.75	1.1	1.1	1.1
整粒电机功率 Granule sizing motor power(kw)		3.0	3.0	3.0	4.0
筛网规格 Screen specification(mm)		d0.6-3	d0.6-3	d0.6-3	d0.6-3
产能 Production capacity(kg/h)		450	450	450	700
重量 Weight(kg)		1000	1100	1200	1600
外形尺寸 Outline dimensions	ΦD	1204	1478	1580	1700
	H	根据工艺要求设计 Designed according to process requirements			

注：本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.

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
Application

The machine is mainly used in the pharmaceutical industry to pellet the granules after granulation, and milling into uniform granules according to the process requirements. It is also widely used in such industries as medicine, chemical industry, food, etc.


Work principle

During operation, the materials are sucked into the buffer tank through the vacuum feeder. The granules flow from the buffer tank to the cavity of the granulator, and the granules in the cavity are subjected to impact, extrusion, shearing, and then discharged through the screen hole, and flow through the discharge port of the guide tube to the hopper. In the closed premise to complete the material transfer, the whole process

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
Features

- ▲The machine height of lifting or falling. Meet clients' requirements.
- ▲The whole process in a closed state to complete,to avoid dust flying,no pollution.
- ▲Meet GMP requirements.

技术参数 Technical parameter

项目 Item	型号 Model	ZLZ-100	ZLZ-200	ZLZ-300	ZLZ-500	ZLZ-700
提升电机功率 Lifting motor power(kw)		1.5	2.2	2.2	3.0	3.0
整粒功率 Granule sizing power(KW)		1.1	3.0	3.0	3.0	4.0
筛网规格 Screen specification(mm)		d0.6-3	d0.6-3	d0.6-3	d0.6-3	d0.6-3
产能 Production capacity(kg/h)		100	450	450	450	700
外形尺寸 Outline dimensions	H	根据用户要求设计 Designed according to user requirements				
	H1					
	L					

注：本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.

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Application

Pharmaceutical: raw material pulverization; wet material milling; dry material milling; unqualified drugs which need to be recycled, can be granulated according to the size of the granules; Food: biscuit crushing and recycling; pulverization of coarse raw materials, the pulverization and granulation of the whole granulated raw material.

Work principle


The raw material is added to the hopper of the pulverizer, and is adjusted to fall into the conical working chamber through the plug-in valve. The rotary granulating knife acts on the raw material to swirl, and the centrifugal force is used to smash the granules toward the screen surface. And at the same time, due to the high speed rotation of the rotary knife and the shearing on the screen surface, and the granules are pulverized into small granules and discharged through the mesh holes. The size of the pulverized granules is adjusted by the mesh of the screen, the distance between the rotary knives and the screen and the speed of the slewing.

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Features

- ▲It thoroughly solves the problem of scraps polluting materials due to the wear of the screen in the oscillating granulator.
- ▲The granule size is easy to control.
- ▲There is no retention dead corner,easy to clean and no overheating of materials.
- ▲Little dust, no vibration and low noise.
- ▲High efficiency, energy saving and high production capacity.
- ▲With beautiful appearance, it is easy to clean.
- ▲With convenient movement, widely application range.

技术参数 Technical parameter

项目 Item	型号 Model	FZM-300	FZM-450	FZM-700	FZM-1000
电机功率 Motor Power(Kw)		2.2	3.0	4.0	5.5
筛网规格 Screen specification(mm)		d0.6-3	d0.6-3	d0.6-3	d0.6-3
产能 Production capacity(kg/h)		300	450	700	1000
外形尺寸 Outline dimensions	H	1300	1450	1650	1875
	H1	897	954	1104	1254
	H2	650	650	700	700
	L	810	920	1000	1100
	W	445	445	500	600

注：本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.

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
Application

Roller compactor adopts continuous feeding and discharging, integrates the extrusion, crushing and granulating functions, directly makes powder into granules. It is especially suitable for granulation of materials which are wet, hot, easily broken down or agglomerated. It has been widely used in pharmaceutical, food, chemical and other industries. In the pharmaceutical industry, granules made by the roller compactor can be directly pressed into tablets or filled into capsules.

Work principle


Mixed with the ingredients of the dry powder material, by the vacuum feeder to the soli material, through the screw conveyor to the pressure chamber, the two high-pressure extrusion wheel to press the material into sheet with the required density and thickness, then through the milling system to make granules meet requirements to complete the granulation process.

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Features

- ▲3 functions in one unit. These are compress.
- ▲Simplified process to obtain highly stable granules. No water or other liquid is added and no need drying, so the loss of active ingredients due to decomposition can be eliminated.
- ▲Pressure rollers with a water cooling system.
- ▲Overall modular design, compact structure, better space utilization.
- ▲The power area is separated from the work area to prevent cross contamination.
- ▲Horizontal feeding structure, exposed installation of extrusion system, easy to clean and maintain.
- ▲Meet GMP requirements.

技术参数 Technical parameter

项目 Item	型号 Model	LG-5	LG-15	LG-50	LG-100	LG-200
送料电机功率 Feeding motor power(kw)		0.37	0.55	0.75	2.2	4
挤压机功率 Extruding motor power(kw)		0.55	0.75	1.5	3	5.5
整粒电机功率 Granulating motor power(kw)		0.37	0.37	0.55	1.1	1.5
油泵电机功率 Oil pump motor power(kw)		0.55	0.55	0.55	0.55	0.55
冷水机功率 Water cooler power(kw)		2.2	2.2	2.2	2.2	2.2
产能 Production capacity(kg/h)		5	15	50	100	200
重量 Weight(kg)		500	700	900	1100	2000

注：本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.

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
Application

The machine is provided with such functions as automatic lifting, mixing, lowering, etc. A bin mixer and variety mixing bin of different specifications ,it can meet the mixing requirements of different batches and diversified varieties of products. It is an ideal machine for mixing in pharmaceutical plants. It is also widely used in such industries as medicine, chemical industry, food, etc.


Work principle

The machine is composed of two machine bases, rotation frame, rotation system, lifting system, control system, etc. Push the mixing bin into the rotation frame and complete the positioning. Start the control system to lift the rotation frame to the mixing height and clamp it reliably. The driving system will carry out automatic mixing according to the set time, rotation speed and other data. After the mixing operation is finished, the rotation frame will stop vertically, the bin will automatically fall to the ground, the machine will stop and process data will be printed. The bin is transferred to the next procedure.

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
Features

- ▲This machine is a new model that our company has extensively absorbed and digested foreign advanced technology, successfully developed in combination with national conditions. Reasonable structure, stable performance, convenient operation, no dead corners, and no exposed bolts.
- ▲The rotating body (mixing bin) symmetry plane is at an angle of 30° to the rotation axis. The materials in the mixing bin rotate along with the rotating body and make tangential movement along the hopper walls, producing strong turning and high-speed tangential movement and achieving the best mixing effect.
- ▲It adopts PLC automatic control . It is provided with the infrared safety device and the discharging butterfly valve with the misoperation preventer to ensure safe production.
- ▲One machine can be equipped with bins of various specifications and quantities.
- ▲The machine effectively controls dust pollution and cross contamination, reduces loss of materials, controls material layering, optimizes production processes.
- ▲Meet GMP requirements.

技术参数 Technical parameter

项目 Item	型号 Model	ZTH-400	ZTH-600	ZTH-800	ZTH-1000	ZTH-1200	ZTH-1500	ZTH-2000
料斗容积 Bin volume(L)		400	600	800	1000	1200	1500	2000
最大装料容积 Maximum loading volume(L)		320	480	640	800	960	1200	1600
最大装料重量 Maximum loading weight(kg)		200	300	400	500	600	750	1000
主轴转速 Spindle rotation speed(rpm)		3-20	3-20	3-20	3-15	3-15	3-15	3-15
回转电机功率 Rotation motor power(kW)		5.5	5.5	5.5	5.5	5.5	7.5	11.0
提升电机功率 Lifting motor power(kW)		1.5	1.5	1.5	2.2	3.0	4	4
整机重量 Machine weight(kg)		1800	2500	2800	3000	3200	3600	4200
外形尺寸 Outline dimensions	H	1780	1780	1880	2070	2150	2240	2410
	H1	1390	1390	1590	1730	1800	1900	2070
	H2	2260	2260	2460	2660	2800	2970	3320
	H3	2710	2710	2940	3200	3350	3500	3810
	L	3290	3290	3290	3660	3710	3910	4010
	W	2150	2150	2150	2300	2300	2300	2300
	W1	1650	1650	1650	1650	1650	1650	1650
	W2	2360	2360	2700	2940	3100	3200	3480

注: 本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.

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Application

The HLT single-column lifting bin mixer is widely used for mixing different components of powders with powders, powders with granules and granules with granules in the production procedure in the pharmaceutical industry. The machine is provided with such functions as automatic lifting, mixing, falling, etc. It can be equipped with various specifications of bin for mixing operation. It is suitable for mixing materials of different varieties and different batches in large and medium-sized pharmaceutical enterprises according to process requirements. Multiple purposes can be realized in one machine.

Work principle

Push the mixing bin into the gyrator arm and lock the nut. Start the control system will carry out automatic mixing according to the set time, rotation speed and other data. After the mixing operation is finished, it is positioned at the horizontal position. The bin will automatically fall to the ground, the machine will stop and process data will be printed. Loosen the locknut of the gyrator arm, and push the bin out and transfer it to the next procedure.

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
Features

- ▲It adopted PLC automatic control; easy control and convenient adjustment of process parameters.
- ▲A single machine can be equipped with bins of different specifications. Improve the efficiency of bin blender, avoid cross contamination
- ▲The mixing barrel can serve both as the dispensing barrel for the previous procedure and as the feeding bin for the next procedure (tablet pressing, filling).
- ▲The mixing uniformity is over 99%, the volume loading coefficient is 0.80 and the barrel is clean, no dead corner .Meet GMP requirements.

技术参数 Technical parameter

项目 Item	型号 Model	HLT-100	HLT-100A	HLT-200	HLT-200A	HLT-300	HLT-300A	HLT-400	HLT-400A	HLT-600	HLT-600A	HLT-800	HLT-800A	HLT-1000	HLT-1000A	HLT-1200	HLT-1200A
料斗容积 Bin volume(L)		100		200		300		400		600		800		1000		1200	
最大装料容积 Maximum loading volume(kg)		80		160		240		320		480		640		800		960	
最大装料重量 Maximum loading weight(kg)		50		100		150		200		300		400		500		600	
主轴转速 Spindle rotation speed(rpm)		3-20		3-20		3-20		3-15		3-15		3-15		3-12		3-10	
整机重量 Machine weight (kg)		1000		1200		1400		1600		1800		2000		2200		2400	
总功率 Total power(kw)		3.7		4.5		4.5		7		8.5		8.5		8.5		8.5	
外形尺寸 Outline dimensions (mm)	L	2550		2830		2870		3000		3330		3380		3430		3500	
	H	2525	2850	2625	2850	2725	2900	2775	2950	2825	3000	2875	3050	2925	3050	3025	3100
	H1	700		700		760		800		960		1060		1060		1160	
	H2	1040		1040		1100		1170		1400		1460		1500		1610	
	H3	1450	1900	1450	1900	1500	1960	1580	2000	1720	2160	1840	2260	1840	2260	1940	2360
	W	740		740		740		740		740		740		740	820	740	820
	W1	1480		1680		1800		1940		1340		2540		2660		2820	

注：本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.

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Application

This machine is developed in combination with the characteristics of a variety of mixers and is suitable for a wide range of mixers . It fully meets GMP requirements and other relevant industry standards. It is especially suitable for the mixing process of products with relatively single varieties and large batch sizes. Such as: pharmaceutical, food, dairy, chemical reagents and other solid preparation industries.

Work principle


The materials enter the mixing barrel through manual feeding, vacuum feeding, floor material dropping or special lifting machine. the symmetrical axis of the hopper is at an angle with the rotation axis of the rotating shaft, The materials of different components will make three-dimensional movement in the closed hopper, producing strong turning, diffusing and shrinking actions and achieving the best mixing effect.

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Features

- ▲With novel design, compact structure and beautiful appearance the machine can reach mixing uniformity 99% and volume loading coefficient is 0.8.
- ▲It has such features as high rotation height, stable running, reliable performance and easy operation.
- ▲The inner and outer walls of the bins are subjected to mirror polishing ,easy discharging and cleaning , no dead corners and without cross contamination.
- ▲Meets GMP requirements.

技术参数 Technical parameter

项目 Item	型号 Model	HF-600	HF-800	HF-1000	HF-1200	HF-1500	HF-2000	HF-3000	HF-4000	HF-6000
容积 Volume(L)		600	800	1000	1200	1500	2000	3000	4000	6000
净负载 Net load(kg)		300	400	500	600	750	1000	1500	1800	3000
转速 Rotation speed(rpm)		2-20	2-16	2-16	2-12	2-12	2-12	2-8	2-8	2-8
电机功率 Motor Power (KW)		3	4.0	4.0	5.5	5.5	7.5	11.0	11.0	15.0
整机重量 Machine weight(kg)		800	850	1200	1650	1650	2000	3000	3500	6000
外形尺寸 Outline dimensions	H	2194	2462	2585	2761	3027	3135	3446	4360	4851
	H1	2160	2300	2400	2345	2546	2576	2848	3813	3971
	H2	1361	1496	1563	1645	1846	1837	2048	2813	2933
	H3	700	800	800	800	900	800	900	1500	800
	L	2608	3075	3188	3479	3618	3869	4120	4470	5223
	L1	590	710	750	780	780	830	800	800	820
	W	1361	1382	1395	1500	1500	1627	1600	2200	2700

注：本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.

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Application

It is mainly used for the tests of mixing different components of dry powders with dry powders, granules with granules, or powders with granules. It is an ideal machine for exploring the best process parameters in material mixing, and an advanced mechanically and electrically integrated experimental machine.

Work principle

When is works,put the materials into the mixing bin. As the symmetrical axis of the bin symmetry plane forms at an angle with the rotation axis, different components of materials will move in the closed hopper, producing strong turning, diffusing and shrinking actions and achieving the best mixing effect.

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
Features

- ▲With novel design, compact structure and beautiful appearance, the machine can reach mixing uniformity 99% and volume loading coefficient is 0.8.
- ▲It has such features as high rotation height, stable running, reliable performance and convenient operation.
- ▲The inner and outer walls of the bin are subject to mirror polishing with easy discharging and cleaning and without any dead corners and cross contamination. Meet GMP requirements.

技术参数 Technical parameter

项目 Item	型号 Model	HLS-5	HLS-10	HLS-15	HLS-20	HLS-30	HLS-50	HLS-100	HLS-150
桶容积 Bin volume(L)		5	10	15	20	30	50	100	150
净负载 Net load(kg)		2.5	5	7.5	10	15	25	50	75
混合转速 Mixing rotation speed(rpm)		3-30	3-30	3-30	3-30	3-25	3-25	3-20	3-20
总功率 Total power(kW)		0.37	0.37	0.37	0.37	0.55	0.55	0.75	1.1
机器重量 Machine weight(kg)		150	160	170	180	190	200	250	280
外形尺寸 Outline dimensions	H	1200	1200	1200	1200	1200	1200	1260	1360
	H1	1030	1030	1030	1030	1030	1030	1110	1210
	H2	650	580	550	520	500	440	440	480
	H3	1050	1120	1150	1180	1200	1280	1430	1600
	L	870	1030	1040	1050	1100	1150	1400	1620
	L1	800	800	800	800	800	800	1000	1050
	W	600	600	600	600	600	600	800	800
	W1	400	600	650	700	730	820	1000	1100

注：本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.

BGB-F 新型高效包衣机+片剂翻斗加料机
High-efficiency coating machine + tablet tipping feeder

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Application

The new type auto coating machine is designed for coating process solutions for all kinds of tablets, pills, granules, which can promote the coating production efficiency and guarantee material stability. Reliable Process, interchangeable pan design, CIP design and good appearance conform to the requirements of GMP. Varieties of functions make the batch production time shorten, satisfying different process requirements. It is the best choice for film coating, sugar coating.

Work principle

Under the rotation action of the coating pan, the tablets or pills to be coated make continuous movement in the pan. The peristaltic pump conveys the coating solution to the spray gun for spraying them on the surfaces of the tablets or pills to be coated. Under negative pressure, the inlet air handling unit supplies clean hot air according to the set procedures and process parameters to dry the tablets or pills. The hot air is discharged from the outlet air handling unit, so that the coating solution sprayed on the tablets or pills rapidly form solid, dense, flat and smooth film, thus the coating process finish.

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Features

- ▲Overall open design on the front side.
- ▲Easy cleaning,no dead corner.
- ▲Wide install area and maintenance space,and pan replacement.
- ▲Unique airbag seal structure,fully enclosed chamber.


新型高效包衣机 High-efficient coating machine

翻斗加料 Tipping feeding

- ▲实现片剂平滑过渡到锅内；
- ▲防止高度差撞击产生碎片、磕片；
- ▲自动提升、斜式翻转，原料可均匀慢速滑入锅内。
- ▲Achieve a smooth transition of the tablet into the pan;
- ▲Prevent fragments and bracts because of the height difference colliding;
- ▲Automatic lifting and tilting, the raw material can slide into the pan evenly and slowly.

喷枪 Spray gun

- ▲独立一体式设计；
- ▲内置管路；
- ▲喷射角度和高度可调。
- ▲Independent and integrated design;
- ▲Built-in pipeline;
- ▲The spray angle and height are adjustable.

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- ▲触控屏工业电脑，多语言操作系统；
- ▲独立模组控制，独立的操作器适合不同要求操作放置；
- ▲在线流量、温度、转速、风速、负压的监控；
- ▲多组程序控制；
- ▲设备和产品信息储存量大，在线打印及输出及上传；
- ▲MES系统端口。
- ▲Touch screen industrial computer, multi-language operating system;
- ▲Independent module control, independent operator is suitable for different requirements to operate;
- ▲Monitoring of online flow, temperature, speed, wind speed and negative pressure;
- ▲Multiple groups of program control;
- ▲Large information storage, online printing and output and upload;
- ▲MES system port.

技术参数 Technical parameter

项目 Item	型号 Model	BGB-75F	BGB-150F	BGB-250F	BGB-350F	BGB-600F
生产能力 Production capacity(kg)		75-50-20	150-75-50	250-150-75	350-250-150	600-350
电机功率 Motor power(kw)		1.5	2.2	3	4	5.5
转速 Rotation speed		2-20	2-15	2-15	2-15	2-10
热风机功率 Hot air fan power(kw)		1.1	1.1	1.5	2.2	5.5
排风机功率 Exhaust fan power(kw)		4	5.5	7.5	11	15
蠕动泵功率 Peristaltic pump power(kw)		0.04(BT100L)	0.05(BT300L)	0.05(BT300L)	0.05(BT300L)	0.05(BT600L)
搅拌桶功率 Mixing barrel power(kw)		2	6	9	9	9
耗蒸汽量 Steam consumption(kg/h)		95	95	225	225	252
耗气量 Air consumption(m³/min)		1.18	1.45	2.03	2.03	2.57

注：本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.

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Application

The BGB-D coating machine is mainly used in medicine and food industries. Meeting GMP requirements, it is a kind of mechanically and electrically integrated coating equipment with high efficiency, energy conservation, good safety and good cleanliness for organic film coating, water soluble coating, dripping pill coating, sugar coating, chocolate and candy coating .

Work principle


Under the rotation action of the coating roller, the tablets or pills to be coated make continuous movement in the roller. The peristaltic pump conveys the coating media to the spray gun for spraying them on the surfaces of the tablets or pills to be coated. Under negative pressure, the inlet air handling unit supplies clean hot air to the tablet bed according to the set procedures and process parameters to dry the tablets or pills. The hot air is discharged from the outlet air handling unit at the bottom of the layer of tablets or pills, so that the coating media sprayed on the tablets or pills rapidly form solid, dense, flat and smooth film to complete the coating operation.

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Features

- ▲The control system consisting of PLC and HMI is designed reasonably with flexible programming, meeting various pharmaceutical process requirements with reliable work and stable performance, conforming to GMP requirements.
- ▲The streamlined diversion plate improve the rate of finished products, reducing broken edges. The narrow upper surface of the diversion plate design improve the adhesion of excipient, saving the excipient and improving medicine quality.
- ▲The constant pressure variable peristaltic pump is used to replace the return flow. The rotation radius of the pan changes with the change of pressure, the output solution and spray volume are automatically balanced, the spray system is optimized, the blocking of the spray gun is avoided, the excipient is saved. Easy to clean without any dead corners.
- ▲It adopts the new type spray gun with even atomizing, large spraying area, adjustable direction. The nozzle is not affected by the loading volume, the anti-blocked device can ensure continuous coating, shorten coating time and save coating excipient.
- ▲Quantitative display and control of technical parameters are realized. SOP operation can be compiled to guide the operation by using quantitative indicators and ensure the consistency of quality of each batch of medicines, fully meeting GMP requirements.

技术参数 Technical parameter

项目 Item	型号 Model	BGB-75D	BGB-150D	BGB-250D	BGB-350D	BGB-600D
生产能力 Production capacity(kg)		75	150	250	350	600
电机功率 Motor power(Kw)		1.5	2.2	3	4	5.5
转速 Rotation speed		2-17	2-15	2-15	2-15	2-6
热风机功率 Hot air fan power(Kw)		1.1	1.1	1.5	2.2	5.5
排风机功率 Exhaust fan power(Kw)		4	5.5	7.5	11	15
蠕动泵功率 Peristaltic pump power(KW)		0.04(BT100L)	0.04(BT300L)	0.05(BT300L)	0.05(BT300L)	0.05(BT600L)
搅拌桶功率 Mixing barrel power(KW)		2	6	9	9	9
电加热功率 Electric heating power(KW)		25	35	39	49	56
耗蒸气量 Steam consumption(Kg/h)		95	95	225	225	252
耗气量 Air consumption(m³/min)		1.18	1.45	2.03	2.03	2.57
主机重量 Weight of main machine(kg)		550	900	1250	1650	2000
外形尺寸 Outline dimensions(L×W×H)(mm)		1200×980×1853	1570×1250×2351	1800×1434×2458	2050×1670×2570	2050×2064×2680
清洗站 Cleaning station	功率 Power(kw)	3	3	3	4	4
	耗水量 Water consumption(T/h)	10	10	10	20	20

注：本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.

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Application

The BGB-C coating machine is mainly used in medicine and food industries. Meeting GMP requirements, it is a kind of mechanically and electrically integrated coating equipment with high efficiency, energy conservation, good safety and good cleanliness for organic film coating, water soluble coating, dripping pill coating, sugar coating, chocolate and candy coating.

Work principle

Under the rotation action of the coating pan, the tablets or pills to be coated make continuous movement in the pan. The peristaltic pump conveys the coating solution to the spray gun for spraying them on the surfaces of the tablets or pills to be coated. Under negative pressure, the inlet air handling unit supplies clean hot air according to the set procedures and process parameters to dry the tablets or pills. The hot air is discharged from the outlet air handling unit, so that the coating solution sprayed on the tablets or pills rapidly form solid, dense, flat and smooth film to complete the coating operation.

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
Features

- ▲The control system consisting of PLC and HMI is designed reasonably with flexible programming, meeting various pharmaceutical process requirements with reliable work and stable performance and conforming to GMP requirements.
- ▲It adopts the new type spray gun with even atomizing, large spraying area, adjustable direction. The nozzle is not affected by the loading volume, the anti-blocked device can ensure continuous coating, shorten coating time and save coating excipient.
- ▲Quantitative display and control of technical parameters are realized. SOP operation can be compiled to guide the operation by using quantitative indicators and ensure the consistency of quality of each batch of medicines, fully meeting GMP requirements.

技术参数 Technical parameter

项目	型号	BGB-75C	BGB-150C	BGB-250C	BGB-350C	BGB-600C
生产能力 Production capacity(kg)		75	150	250	350	600
电机功率 Motor Power(Kw)		1.5	2.2	3	4	5.5
转速 Rotation speed		2-17	2-15	2-15	2-15	2-6
热风机功率 Hot air fan power(Kw)		1.1	1.1	1.5	2.2	5.5
排风机功率 Exhaust fan power(Kw)		4	5.5	7.5	11	15
蠕动泵功率 Peristaltic pump power(KW)		0.04(BT100L)	0.04(BT300L)	0.05(BT300L)	0.05(BT300L)	0.05(BT600L)
搅拌桶功率 Mixing barrel power(KW)		2	6	9	9	9
电加热功率 Electric heating power(KW)		25	35	39	49	56
耗蒸汽量 Steam consumption(Kg/h)		95	95	225	225	252
耗气量 Air consumption(m³/min)		1.18	1.45	2.03	2.03	2.57
主机重量 Weight of main machine(kg)		550	900	1250	1650	2000
外形尺寸 Outline dimensions(L×W×H)(mm)		1200×980×1853	1570×1250×2351	1800×1434×2458	2050×1670×2570	2050×2064×2680

注：本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.

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Application

The BGB-F coating machine is mainly used in medicine and food industries. Meeting GMP requirements, it is a kind of mechanically and electrically integrated coating equipment with high efficiency, energy conservation, good safety and good cleanliness for organic film coating, water soluble coating, dripping pill coating, sugar coating, chocolate and candy coating.

Work principle

Under the rotation action of the coating roller, the tablets or pills to be coated make continuous movement in the roller. The peristaltic pump conveys the coating media to the spray gun for spraying them on the surfaces of the tablets or pills to be coated. Under negative pressure, the inlet air handling unit supplies clean hot air to the tablet bed according to the set procedures and process parameters to dry the tablets or pills. The hot air is discharged from the outlet air handing unit at the bottom of the layer of tablets or pills, so that the coating media sprayed on the tablets or pills rapidly form solid, dense, flat and smooth film to complete the coating operation.

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
Features

- ▲Changeable coating pans can be replaced according to the different batch size.
- ▲Easy operation to replace the coating pan.
- ▲With compact and intergrated structure, air intake at the clean area.

技术参数 Technical parameter

项目 Item	型号 Model	BGB-5F	BGB-10F	BGB-20F	BGB-40F
生产能力 Production capacity(kg)		1-3-5	10-5-3	5-10-20	20-40
电机功率 Motor Power(Kw)		0.37	0.55	0.55	1.10
转速 Rotation speed		2-30	2-30	2-26	2-26
热风机功率 Hot air fan power(Kw)		0.4	0.4	0.4	0.75
排风机功率 Exhaust fan power(Kw)		0.75	1.5	2.2	3
蠕动泵功率 Peristaltic pump power(KW)		0.04(BT100L)	0.04(BT100L)	0.04(BT100L)	0.04(BT100L)
搅拌桶功率 Mixing barrel power(KW)		/	1	1	2
电加热功率 Electric heating power(KW)		5	10.5	10.5	14.7
耗蒸气量 Steam consumption(Kg/h)		67	90	90	90
耗气量 Air consumption(m³/min)		0.45	0.65	0.65	0.65
主机重量 Weight of main machine(kg)		450	450	500	600

注：本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.

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Application

The containment coating machine is mainly suitable for tablet coating of API materials. The isolation system ensures the safety of operators, materials and environment, and realizes multiple coating process.

Work principle

Under the rotation action of the coating pan, the tablets or pills to be coated make continuous movement in the pan. The peristaltic pump conveys the coating solution to the spray gun for spraying them on the surfaces of the tablets or pills to be coated. Under negative pressure, the inlet air handling unit supplies clean hot air according to the set procedures and process parameters to dry the tablets or pills. The hot air is discharged from the outlet air handling unit, so that the coating solution sprayed on the tablets or pills rapidly form solid, dense, flat and smooth film, thus the coating operation finish.

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

Features

- ▲PLC programmable controller, automatic control program;
- ▲HMI interface, easy operation;
- ▲Three-level administrator password;
- ▲Remote communication function (optional);
- ▲Data export and print function;
- ▲Automatic operation& manual operation mode;
- ▲Temperature and humidity and differential pressure monitoring function;
- ▲Pressure alarm function (sound and light signal);
- ▲Lighting: illuminance ≥ 3001x;
- ▲Reserve detection port of dust particle, floating bacteria and sedimentation bacteria;
- ▲Reserve port for sterilization;
- ▲CIP cleaning;
- ▲Meet the OEB4 standard;
- ▲Push-push filter for getting rid of dust.

技术参数 Technical parameter

项目 Item	型号 Model	BGB-5M	BGB-10M	BGB-20M	BGB-40M
生产能力 Production capacity(kg)		1-3-5	10-5-3	5-10-20	20-40
电机功率 Motor Power(Kw)		0.37	0.55	0.55	1.10
转速 Rotation speed		2-30	2-30	2-26	2-26
热风机功率 Hot air fan power(Kw)		0.4	0.4	0.4	0.75
排风机功率 Exhaust fan power(Kw)		0.75	1.5	2.2	3
蠕动泵功率 Peristaltic pump power(KW)		0.04(BT100L)	0.04(BT100L)	0.04(BT100L)	0.04(BT100L)
搅拌桶功率 Mixing barrel power(KW)		/	1	1	2
电加热功率 Electric heating power(KW)		5	10.5	10.5	14.7
耗蒸气量 Steam consumption(Kg/h)		67	90	90	90
耗气量 Air consumption(m³/min)		0.45	0.65	0.65	0.65
主机重量 Weight of main machine(kg)		450	450	500	600

注：本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.

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Application



It is mainly used in the processing of pelletized varieties that require layer coating. At present, various types of regular or irregular coating processing with core diameters of 0.02~10mm can be processed, and the core layer is wrapped with adhesive. The realization process of the coating layer can realize a one-time input of a certain amount of core products and realize a fully automatic coating process. At present, there are more applications in the agricultural seeds, food, and chemical industries. The quality of automatic coating pill products should not only meet the quality standards, but also further improve the appearance, strength, cracking and other indicators. Automatic coating of pelletizing equipment to achieve seed formation.

Work principle

Starting the automatic coating process after adjusting the spray gun angle and flow amount. The core is subjected to high-speed centrifugation to form an outer ring circulation. Under the action of the baffle, the core is flipped and dropped to the middle ring and the inner ring. The adhesive is sprayed on the falling parabola, combined with the added coating powder, quickly wrapped. During the centrifuged process, the core and the coating powder and the binder are subjected to the centrifugal force to compact and form firm coating. The required finished particle size is achieved after numerous cycles.

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Features

- ▲ Automatically complete the expected pelleting size by one-time feed;
- ▲ Ensure core integrity, uniformity, and individual independence;
- ▲ Save core and excipient, low production costs;
- ▲ Improve the cracking and performance;
- ▲ For the mini and irregular core, the overall volume can be increased and the shape and size are uniform after the pelletization treatment, which is good to mechanized operation.

技术参数 Technical parameter

项目 Item	型号 Model	WH-5	WH-10	WH-20	WH-50
产能(核心)批次 Capacity (core) batch(kg/批batch)		0.5	1	2	5
核心规格 Core specification(mm)		0.02~10			
增量倍数 Incremental multiple		< 100		< 60	
容量 Volume(L)		240	300	400	600
喷枪 Spray gun(PC)		2	2	2	3
功率 Power(kw)		5	7	9	15
外形尺寸 Outline dimensions(mm)		900 × 600 × 2000	1100 × 700 × 2500	1200 × 800 × 2700	1500 × 900 × 2700

注：本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.

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Application

This machine is mainly used for conveying, lifting and feeding of solid materials such as powder, granule and flake in pharmaceutical, food and chemical industries. It is supporting equipment for blenders, tablet presses, coating machines, capsule filling machines, packaging machines, etc.

Work principle


The equipment mainly consists of a fixed base, a lifting system, a rotating system and a shaft-turning system. The servo motor is connected to the geared motor as the driving source, and the motor brake is used to achieve high-speed safe movement and positioning, and it is equipped with anti-falling devices. One-touch operation can be realized by setting the movement route according to the program.

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Features

- ▲Automatic completion of clamping-lifting-rotating-down-discharging-resetting process by one-touch button operation.
- ▲Driving power is driven by stepper motor instead of hydraulic oil;
- ▲The motor with brake protects in case of power failure;
- ▲The air expansion shaft is used to fix the barrel and the lifting fork arm;
- ▲start and stop freely at any position by limit position device;
- ▲Stainless steel belt is used to replace the column curtain;
- ▲Linear guide rail, more resistance to torque, smoother linear motion;
- ▲The double-row chain structure improves the load and makes the operation more stable;
- ▲Using the medium-sealed solenoid valve, real-time pressure detection, it plays the role of accidental gas-breaking protection, and the gas-up axis can guarantee the loading state;
- ▲add the fork arm rotation structure, precise docking, also be used as a single-column blender (fork-arm rotary lifting blender).

技术参数 Technical parameter

型号 Model	净负载 Net load(kg)	电机功率(KW) Motor power	驱动轴 Drive axis		外形尺寸 Outline dimensions (mm)
ATG-200	200	7	2轴(标配) 提升、主臂旋转 2-axis (standard) lift, main arm rotation	3轴(选配) 提升、主臂旋转、载臂旋转 3-axis (optional) lift, main arm rotation, arm rotation	根据用户要求设计 Designed according to user requirements
ATG-400	400	12			
ATG-600	600	18			
ATG-800	800	22			

注：本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.

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Application

The machine is mainly used for conveying and charging solid materials in the pharmaceutical industry. It can work with the mixer, tablet press, capsule filler, etc. It is also widely used in such industries as medicine, chemical industry, food, etc.

Work principle


The machine is mainly composed of chassis, column, lifting system, etc. When it works, push the hopper loaded with materials into the lifting crotch of the hopper lifter. Press the lifting button and the hopper will make lifting movement. After the hopper reaches the designated position, rotate the chassis to realize closed connection with the charging equipment. Start the discharging butterfly valve to transfer the materials to the next procedure.

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Features

The machine is a new machine researched and developed successfully by our company according to China's actual conditions after absorbing and digesting internationally advanced technology. It has such features as reasonable structure, stable performance, convenient operation, no dead corners, and no exposed bolts. The machine is easy to clean, effectively controls dust pollution and cross contamination, optimizes production processes, and fully meets GMP requirements for medicine production.

技术参数 Technical parameter

型号 Model	净负载 Net load(kg)	液压站容积 Hydraulic station volume	液压站功率 Hydraulic station power	H	H1	L
YTG-200	200	21L	2.2kw	根据用户要求设计 Designed according to user requirements		
YTG-300	300					
YTG-400	400					
YTG-500	500					
YTG-600	600	31L or 54L	3.0kw			
YTG-800	800					
YTG-1000	1000					
YTG-1200	1200					
YTG-1500	1500					

注：本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.

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Application


The machine is mainly used for conveying and charging solid materials in the pharmaceutical industry. It can work with the mixer, granule mill machine, tablet press, coating machine, capsule filling machine, etc. It is also widely used in medicine, chemical, food industries and so on.

Work principle

The machine is mainly composed of chassis, column, lifting system, etc. When it works, push the bin loaded with materials into the lifting fork of the lifter, start the lifting button and the bin will make lifting movement. After the bin reaching the designated position, rotate the chassis to realize closed connection with the charging equipment. Start the discharging butterfly valve to transfer the materials to the next process.

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Features

The PLC program is used to control the servo motor drive, which has strong anti-overload capability, stable operation at low speed, good controllability, fast response, high sensitivity, and significantly reduces heat and noise. The new type equipment is concise and easy to operate. The bin discharging butterfly valve is easy to assemble, disassemble and clean, meeting GMP requirements in the pharmaceutical industry. The equipment has the advantages of light weight, large load and good running stability. The cart wheels are large with low rolling resistance and easy movement. It has multiple purposes, saving time, saving labor and reducing environmental pollution. It has low maintenance requirements, which can reduce the maintenance cost.

技术参数 Technical parameter

型号 Model	净负载 Net load(kg)	液压站容积(L) Hydraulic station volume	液压站功率(kw) Hydraulic station power	H	L
YTY-100	100	18	1.5	1950	根据用户要求设计 Designed according to user requirements
YTY-150	150				
型号 Model	净负载 Net load(kg)	电动功率 Electric power(kw)			
YTD-150	150	1.5			

注：本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.

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Application


The machine is mainly used for upward and downward conveying of solid materials (bagged, barreled) at two different clean areas among floors in the pharmaceutical industry. The machine effectively reduces cross contamination of dust, lowers labor intensity and fully meets GMP requirements for medicine production. It is an ideal material conveying machine. It is also widely used in chemical, food, electronics and so on industries.

Work principle

The lifting machine among floors is mainly composed of column, lifting platform system, control system, etc. When it works, first send the materials to the lifting platform, close the safety fence and start the lifting operation lever on the control panel to lift the materials to the designated position, lock the upper limit device, when the lifting platform stops movement, open the safety fence and push out the materials to complete upward and downward conveying of materials among floors.

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Features

- ▲The column structure of the machine is made of solid steel frame. The two pieces of special heavy forklift channel steel have high strength after welding. The exterior surfaces are made of matte stainless steel pattern plate, which is easy to clean.
- ▲The lifting groove is separated by curtain, which meets the requirement of clean area.
- ▲It is equipped with safety electric interlocking and upper and lower limit devices among floors to ensure equipment safety during the movement, fully meeting GMP requirements for medicine production.

技术参数 Technical parameter

型号 Model	净负载 Net load(kg)	液压站容积 Hydraulic station volume	液压站功率 Hydraulic station power	H	H1	W	L
CT-200	200	34L	2.2kw	根据用户要求设计 Designed according to user requirements			
CT-300	300						
CT-400	400						
CT-500	500						
CT-600	600	31L or 54L	3.0kw				
CT-800	800						
CT-1000	1000						
CT-1200	1200						
CT-1500	1500						

注: 本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.

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Application

The ZLXHD series bin cleaning machine is mainly used for cleaning stainless steel drums, IBC transfer bin and mixing transfer bins in medicine, chemical, food and so on industries. The ZLXHD series bin cleaner can effectively clean foreign matters remaining on the inner and outer surfaces of the bin to avoid cross contamination of various ingredients during the production. It is an indispensable machine in pharmaceutical enterprises. It is also one of the necessary equipment for pharmaceutical enterprises to meet the requirements of GMP in the production process of solid pharmaceutical preparations.

Work principle


Push the bin into the cleaning room and position it, close the door and start the cleaner. The cleaner will carry out cleaning, drying and cooling of the bin according to the set procedures. After all procedures are completed, open the cleaning machine door and push out the bin and send it to the storage station for future use.

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Features

- ▲It provides uniform cleaning standard for the cleaning of vessels and easily makes the traceability and certification of the cleaning process.
- ▲It improves production efficiency.
- ▲It reduces labor intensity of workers.
- ▲The equipment is a machine integrating cleaning, drying and cooling functions. It adopts PLC control during the whole course, which is with high degree of automation and simple operation.
- ▲The running parameters are recorded automatically to verify the running results.

技术参数 Technical parameter

项目 Item	型号 Model	ZLXHD-600	ZLXHD-800	ZLXHD-1000	ZLXHD-1200	ZLXHD-1500	ZLXHD-2000
总功率 Total power(Kw)		12.25	12.25	12.25	12.25	12.25	12.25
水泵功率 Pump power(Kw)		4	4	4	4	4	4
水泵流量 Pump flow(t/h)		20	20	20	20	20	20
泵压力 Pump pressure(Mpa)		0.35	0.35	0.35	0.35	0.35	0.35
引风机功率 Induced draft fan power(Kw)		2.2	2.2	2.2	2.2	2.2	2.2
排风机功率 Exhaust fan power(Kw)		5.5	5.5	5.5	5.5	5.5	5.5
蒸汽压力 Steam pressure(Mpa)		0.4-0.6	0.4-0.6	0.4-0.6	0.4-0.6	0.4-0.6	0.4-0.6
蒸汽流量 Steam flow(kg/h)		1300	1300	1300	1300	1300	1300
压缩空气压力 Compressed air pressure(Mpa)		0.4-0.6	0.4-0.6	0.4-0.6	0.4-0.6	0.4-0.6	0.4-0.6
压缩空气耗量 Compressed air consumption(m³/min)		3	3	3	3	3	3
整机重量 Equipment weight(t)		4	4	4.2	4.2	4.5	4.5
外形尺寸 Outline dimensions (mm)	L	2000	2000	2200	2200	2200	2200
	H	2820	3000	3100	3240	3390	3730
	H1	1600	1770	1800	1950	2100	2445
	H2	700	700	700	700	700	700

注：本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.

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Double chamber bin cleaning machine

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
Application

The ZLXHS series double chamber bin cleaning machine is mainly used for cleaning stainless steel drums, IBC transfer bins and mixing transfer bins in medicine, chemical, food and so on industries.
The ZLXHS series double chamber bin cleaning machine can effectively clean foreign matters remaining on the inner and outer surfaces of the bin to avoid cross contamination of various ingredients during the production. It is an indispensable machine in pharmaceutical enterprises. It is also one of the necessary equipment for pharmaceutical enterprises to meet the requirements of GMP in the production process of solid pharmaceutical preparations.

Work principle

The machine is composed of conveying system, cleaning chamber, drying chamber, pumping station, air handling unit, electric cabinet, operation cabinet, etc. The conveying, cleaning, drying and cooling of bin are automatically completed according to the set system. When working, the bin is manually pushed into the conveying system and the bin is automatically conveyed to complete cleaning, drying. The cleaning is completed in the first chamber, and then the automatic conveying system is transferred to the second chamber for drying, the bin Cooling and finally stopping outside the drying room. It is manually removed and stored at the storage station for future production use.

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

Features

- ▲ It provides uniform cleaning standard for the cleaning of vessels and easily makes the traceability and certification of the cleaning process.
- ▲ It improves production efficiency.
- ▲ It reduces labor intensity of workers.
- ▲ The equipment is a machine integrating cleaning, drying and cooling functions. It adopts PLC control during the whole course, which is with high degree of automation and simple operation.
- ▲ The running parameters are recorded automatically to verify the running results.
- ▲ It fully meets GMP requirements.

技术参数 Technical parameter

项目 Item	型号 Model	ZLXHS-600	ZLXHS-800	ZLXHS-1000	ZLXHS-1200	ZLXHS-1500	ZLXHS-2000
总功率 Total power(Kw)		14.25	14.25	14.25	14.25	14.25	14.25
水泵功率 Pump power(Kw)		4	4	4	4	4	4
水泵流量 Pump flow(t/h)		20	20	20	20	20	20
泵压力 Pump pressure(Mpa)		0.35	0.35	0.35	0.35	0.35	0.35
引风机功率 Induced draft fan power(Kw)		2.2	2.2	2.2	2.2	2.2	2.2
排风机功率 Exhaust fan power(Kw)		5.5	5.5	5.5	5.5	5.5	5.5
蒸汽压力 Steam pressure(Mpa)		0.4-0.6	0.4-0.6	0.4-0.6	0.4-0.6	0.4-0.6	0.4-0.6
蒸汽流量 Steam flow(kg/h)		1300	1300	1300	1300	1300	1300
压缩空气压力 Compressed air pressure(Mpa)		0.4-0.6	0.4-0.6	0.4-0.6	0.4-0.6	0.4-0.6	0.4-0.6
压缩空气耗量 Compressed air consumption(m³/min)		3	3	3	3	3	3
整机重量 Equipment weight(t)		6.6	6.6	7	7	7.2	7.2
外形尺寸 Outline dimensions (mm)	L	7000	7000	8100	8100	8100	8100
	H	2820	3000	3000	3240	3390	3730
	W	4100	4100	4100	4100	4600	4600
	H1	1600	1770	1800	1950	2100	2445
	H2	700	700	700	700	700	700

注：本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.


QXY 移动清洗机
Moveable cleaning machine

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Application

With compact structure, convenient operation, good sealing performance and high adaptability, the moveable cleaning machine produced by our company is used for cleaning industrial equipment and other parts that require cleaning in medicine, food and fine chemical industries. It is a processing machine designed and manufactured in accordance with GMP requirements.

Work principle


The pressure in the boosting pump is used to convey the mixture of the cleaning liquid and the desired water source. Through different water source connections, different water inlet valves can be operated according to the needs to connect with water sources. The detergent adding volume is controlled by valve then enters to boosting pump after mixing. Under the action of the boosting pump, flow output is formed within the pressure range of the pump according to the parameters in the pump height-flow performance table. The output flow varies with the pressure.

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Features

- ▲The machine can be equipped with three water pipe inlets for users to select. The water outlets can be connected to tees or crosses to clean several machines simultaneously, which is convenient and easy to operate.
- ▲The machine is equipped with a cleaning nozzle that can turn for 360° and a high-pressure water gun. The procedure of bin cleaning can be completed through the operation by just one person.
- ▲The machine is installed with castors at its bottom, which is easy to move.
- ▲It is featured by concise structure and easy installation and disassembly.

技术参数 Technical parameter

项目 Item	型号 Model	QXY-10	QXY-20
流量 Flow(T/h)		10	20
电机功率 Motor power(kw)		3.0	4.0
压力 Pressure(MPa)		0.62	0.35
设备重量 Equipment weight(kg)		100	125
外形尺寸 Outline dimensions (mm)	L	1490	1490
	W	400	400
	H	1150	1150

注：本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.

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Application

The machine is mainly used for online cleaning of equipment in the production process of solid preparations, used together with equipment such as fluid bed machine, coating machine and wet mixing granulator. It is also suitable for the cleaning of industrial equipment such as pharmaceutical, food, fine chemicals, etc. It is a cleaning equipment designed and manufactured in accordance with GMP requirements.

Work principle

The machine is composed of booster pump, heat exchanger, water storage tank, control system and so on, and the various water source are pressurized by a booster pump to clean the inner surface of the device, Hot water source is generated by tap water heated by the heat exchanger. The cleaning agent is controlled by the valve to add the amount, entering to the booster pump after mixing. The cleaning procedure is set according to the process requirements. After the cleaning is completed, the compressed air is blown according to the program to dry the water in the pipeline, and the cleaned equipment is dried for future production.

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
Features

- ▲This machine is optionally equipped with heat exchanger, if users have no hot water source.
- ▲The outlet can be connected to several devices for online cleaning.
- ▲Cleaning procedures, processes, etc. are set on the touch screen, the automation program is high, and the operation is simple.


技术参数 Technical parameter

项目 Item	型号 Model	ZQX-10	ZQX-20
水泵流量 Pump flow(t/h)		10	20
水泵功率 Pump power(Kw)		3	4
泵压力 Pump pressure(Mpa)		0.62	0.35
蒸汽压力 Steam pressure(Mpa)		0.4-0.6	0.4-0.6
蒸汽流量 Steam flow(kg/h)		1300	1300
设备重量 Equipment weight(kg)		450	475
外形尺寸 Outline dimensions (mm)	L	2400	2500
	W	1075	1090
	H	2100	2100

注：本公司可按用户要求定制生产
Note: Our company can customize products according to user requirements.

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

HZT全自动混合周转料斗
AUTOMATIC MIXING TRANSFER BIN

技术参数 Technical parameter

型号 Model	容积(L) Volume	W	H	H1	H2	t	Φ1	Φ2	DN	重量(kg) Weight
HZT-400	400	800	1270	1210	180	2.5	400	125	200	150
HZT-600	600	1000	1280	1210	130	2.5	400	125	200	185
HZT-800	800	1000	1480	1420	130	2.5	400	125	200	200
HZT-1000	1000	1180	1490	1430	130	3	500	150	200	245
HZT-1200	1200	1180	1640	1580	130	3	500	150	200	260
HZT-1500	1500	1180	1850	1780	130	3	500	150	200	280
HZT-2000	2000	1180	2210	2130	130	3	500	150	200	310


ZT方形周转料桶
SQUARE TRANSFER BIN

技术参数 Technical parameter

型号 Model	容积(L) Volume	W	H	H1	t	Φ1	Φ2	DN	重量(kg) Weight
ZT-200	200	600	1140	185	2	400	125	150	145
ZT-300	300	800	1180	210	2.5	400	125	200	160
ZT-400	400	800	1320	220	2.5	400	125	200	180
ZT-600	600	1000	1450	230	2.5	500	125	200	200
ZT-800	800	1000	1650	230	2.5	500	125	200	220
ZT-1000	1000	1200	1590	200	2.5	500	150	250	262
ZT-1200	1200	1200	1720	200	2.5	500	150	250	280
ZT-1500	1500	1200	1950	200	2.5	500	150	250	315
ZT-1800	1800	1200	2160	200	3	500	150	250	335
ZT-2000	2000	1200	2290	200	3	500	150	250	370

ZTY圆形周转料桶
ROUND TRANSFER BIN

技术参数 Technical parameter

型号 Model	容积(L) Volume	H	t	Φ1	Φ2	Φ3	DN	重量(kg) Weight
ZTY-200	200	1050	2	750	125	400	200	89
ZTY-400	400	1287	2.5	850	125	400	200	126
ZTY-600	600	1457	2.5	1000	125	500	200	178
ZTY-800	800	1530	2.5	1100	125	500	200	190
ZTY-1000	1000	1615	2.5	1250	150	500	250	262

配件
Accessories

KGT料桶卡箍
DRUM CLAMP

料桶用卡箍
Clamp for drum


技术参数 Technical parameter

型号 Model	DN
KGT-500	500

KGF蝶阀卡箍
BUTTERFLY VALVE CLAMP

蝶阀用卡箍
Clamp for butterfly valve

技术参数 Technical parameter

型号 Model	DN
KGF-150	150
KGF-200	200
KGF-250	250

KGD料斗卡箍
BIN CLAMP

料斗用卡箍
Clamp for bin


技术参数 Technical parameter

型号 Model	DN
KGD-400	400
KGD-500	500

KGG 管道卡箍
PIPE CLAMP

管道用卡箍
For pipeline connection

技术参数 Technical parameter

型号 Model	D
KGG-4	119
KGG-6	167
KGG-8	218
KGG-10	269
KGG-12	319

配件
Accessories

LT 不锈钢料桶
STAINLESS STEEL DRUM

配提升加料机及物料周转
For pharma lifter and material transferring

技术参数 Technical parameter

型号 Model	主要结构尺寸 Main structural dimensions	
	Φ D	h
LT-50	400	400
LT-100	500	500
LT-150	585	600

ZCQ 真空出料器
VACUUM DISCHARGER

料斗，容器真空出料
For bin/container vacuum discharging


技术参数 Technical parameter

型号 Model	主要结构尺寸 Main structural dimensions		
	Φ D	Φ D1	Φ D2
ZCQ-150-38	150	38	230
ZCQ-150-50	150	50	230
ZCQ-200-38	200	38	280
ZCQ-200-50	200	50	280
ZCQ-200-63	200	63	280
ZCQ-250-50	250	50	330
ZCQ-250-63	250	63	330

LTC 十字小车
CROSS CART

料桶搬运
For drum conveying

技术参数 Technical parameter

型号 Model	主要结构尺寸 Main structural dimensions		
	Φ D	H	Φ d
LTC-50	500	160	75
LTC-100	600	160	75

YXG Y型管
Y TYPE PIPE

配浆、制浆、保温、压片机等加料
Charging for paste preparing & pulping,
heat insulation tank, tablet press

技术参数 Technical parameter

型号 Model	主要结构尺寸 Main structural dimensions				
	Φ D1	Φ D2	Φ D3	L	H
YXG-150	150	219	100	829	510
YXG-200	200	269	150	829	645
YXG-250	250	269	200	795	480

DJL 硅胶接管
SILICONE CONNECTION PIPE


联接料斗加料
Connecting with bin to charge material

技术参数 Technical parameter

型号 Model	主要结构尺寸 Main structural dimensions		
	Φ D	Φ D1	H
DJL-150	175	145	170
DJL-150	225	195	170
DJL-250	275	245	170

α-β 阀门
VALVE


技术参数 Technical parameter

接口尺寸 Interface size		2.0"、2.5"、3.0"、4.0"、6.0"、8.0"
接口标准 Interface standard		Tri-clamp、PN6/PN10法兰
材质 Material	阀体材质 Valve material	316L(与介质接触部件Material contact parts) 304(不与介质接触部件Remaining parts)
	密封材质 Sealing material	白色氟化橡胶(标配)，符合FDA要求 White fluorinated rubber(Standard), Meet FDA requirements 三元乙丙橡胶、硅橡胶、全氟橡胶(选配) Ethylene-Propylene-Diene Monomer, Silicone rubber, Perfluoro rubber(Optional)
操作方式 Operating mode		手动、自动 Manual, Automatic
密封等级 Seal grade		OEB 4级 Level 4(OEL ≤ 10 μg/m³)
压力等级 Pressure rating		-0.1~+0.5Mpa
灭菌方式 Sterilization method		SIP
防爆等级 Anti-explosion level		ATEX II2 GD T4
选配附件 Optional accessories		主动压力塞、主动保护塞、主动清洗部件、被动压力盖、被动保护盖、被清洗部件 Active pressure plug, Active protection plug, Active cleaning unit, Passive pressure cap, Passive protective cover,leaned part
表面处理 Surface treatment		Ra < 0.4, 标配Standard(与介质接触部件Material contact parts) Ra < 0.8(不与介质接触部件Remaining parts)

DF 蝶阀
BUTTERFLY VALVE

输送，控制，配料斗
Used with bin to discharge and control material

技术参数 Technical parameter

型号 Model	主要结构尺寸 Main structural dimensions		
	Φ D1	Φ D2	H
DF-150	150	230	55
DF-200	200	280	55
DF-250	250	330	60

FLG 分料盖
SUB-FEEDING BIN COVER

配料斗分料
Used with bin to sub-feed material

技术参数 Technical parameter

型号 Model	主要结构尺寸 Main structural dimensions		
	Φ D	Φ D1	Φ d
FLG-400-150	150	400	38
FLG-400-200	200	400	38
FLG-400-250	250	400	38
FLG-500-150	150	500	38
FLG-500-200	200	500	38
FLG-500-250	250	500	38