



HASELMEIER™

A medmix Brand

PiccoJect™

PRODUCT BROCHURE



EXCELLENCE THROUGH
SIMPLICITY

Copyright © 2022 medmix Ltd. All rights reserved.



DRUG DELIVERY





Key features of PiccoJect™

PiccoJect™ is a highly compact, customizable and fully featured two-step autoinjector designed for subcutaneous delivery of drug products.

8


Extremely low part count

Only 8 parts in total – reduces manufacturing and scale up challenges.



Ergonomic shape & small size

Thanks to its compact flat form the PiccoJect™ is easy to grip and hold.




Large wrap-around drug window with customizable size

Facilitates visual inspection of the drug prior to use.



Colored status indicator

Provides easy-to-understand binary information about the usage status of the autoinjector.




Optimized spring forces

A parallel layout allows for larger diameter springs, which enables optimized spring forces.



Audible clicks at the start and end of injection

Help to ensure that the user holds the device in place until the full dose has been injected.



Future proof

Multiple connectivity options to integrate into a digital ecosystem.



Any standard 1 ml long or 2.25 ml pre-filled syringe

Design accommodates glass or plastic syringes with a small round or cut flange.



Committed to sustainability

Low carbon intensity as a standard offering: Utilizing materials with sustainable feedstocks, investment in green electricity, and development of regional supply chains for the US and European market.



A full-service platform

Combination product development, design verification, final assembly, secondary packaging, labeling, and serialization included.



PiccoJect™ – 100

PiccoJect™ – 225

Technical specifications

SPECIFICATIONS	PICCOJECT™ – 100	PICCOJECT™ – 225
Part count	5 plastic components, 2 springs, 1 metal component	
Primary container	1 ml long glass or plastic syringe	2.25 ml glass or plastic syringe
Syringe flange	Small round or cut flange	Small round or cut flange
Fill volume (how to adapt the fill volumes)	0.2 – 1 ml	0.6 – 2 ml
Injection time*	< 10 s	< 15 s
Viscosity**	Up to 20 cP	
Needle insertion depth	6 mm (nominal); customization possible	
User feedback	Audible click at start and end of dose; visual feedback in drug window and dedicated status indicator	
Needle safety	Automatic needle shielding with needle hidden before, during and after use	
Needle type and gauge	27G and 29G, normal wall through special thin wall	
Weight	26 g (without syringe)	34 g (without syringe)
Dimensions	H = 130 (cap on), W = 26, D = 15 mm	H = 134 (cap on), W = 30, D = 19 mm

*, ** Injection time and viscosity capability are dependent on needle diameter and fill volume.



Simplicity is key to efficiency

'Excellence through simplicity' sums up the key features of PiccoJect™.

It all starts with an extremely low part count:
The PiccoJect™ autoinjector is made up of only eight parts.
Apart from the syringe, the highly integrated delivery mechanism consists of five injection-molded plastic parts, two springs and one metal component.



The delivery mechanism of the PiccoJect™ autoinjector accommodates any standard 1 ml long and 2.25 ml pre-filled glass or plastic syringe with a small round or cut flange. The same mechanism is available with two different cross-sections, tailored to the applicable syringe size to provide two discrete and user-friendly device form factors.

only 8 parts:
reduced
manufacturing &
scale up
challenges

ergonomic
shape & large
drug window:
excellent user
experiences

materials with
sustainable
feedstocks &
regional supply
chains




Easy to customize

Thanks to its versatile design, the PiccoJect™ autoinjector can easily adapt to a range of customer requirements. Standard customization options include cap color, the size of the drug window, and spring force. Agility is key: these customization options are available without any detrimental impact on development timelines.



fully customizable label



different color options
for the cap



customizable size of
the drug window



PiccoJect™ is designed for the future

The design includes a number of connectivity options that are either integrated or provided as an add-on module, without modifying the core mechanism or compromising handling.

For connectivity to an app, an add-on bluetooth module is being developed that can be re-used on multiple pens. The module is easy to attach and has an expected two years battery life without the need to recharge.

As the coverage of IoT networks increases and the costs of hardware go down, we can also provide a version of PiccoJect™ that connects directly to the cloud, allowing data to be collected without requiring the use of a patient app.

Automatic collection and storage of treatment data helps patients save time because there is no need to document each injection event on paper.





Committed to sustainability

At Haselmeier, sustainability is embedded in our daily business and in the foundation of our corporate strategy. Our sustainability policy enables us to reduce waste and minimize environmental impact. We keep our supply chain short and simple for reliable deliveries. Excellence through simplicity applies not only to our products and processes but also to the way we collaborate with our customers.

Eco-design principles

The underlying design concept of the PiccoJect™ autoinjector was driven by sustainability requirements. PiccoJect™ is Haselmeier's first product which was developed according to our eco-design principles. This resulted in a design that features low carbon intensity materials as standard for all plastic components.


In line with the medmix approach to keep shipping routes short, our manufacturing and distribution facilities are located in the geographies where our customers and partners are based. The PiccoJect™ autoinjector is manufactured at Haselmeier sites that rely on low-carbon electricity.




Delivering a better user experience

To help ensure adherence, the PiccoJect™ autoinjector is designed from the ground up for ease of use.


- 1 Remove cap**
by pulling it
away from the
autoinjector body




- 2 Align with body** (either
abdomen or thigh)




- 3 Gently press down**
so that the needle
guard depresses
(you may hear or
feel a click)




- 4 Hold still and wait**
for the injection to
finish (drug window
and status indicator
change to yellow)



- 5 Wait a further 3 seconds**
after the indicator
changes to allow the
injection to complete



- 6 Remove the**
autoinjector from
the skin and place
in the sharps bin



Users are excited

Results from a recent study

Scope and setup

- PiccoJect™ was compared in a preference study with comparable marketed devices.
- Users (n=24) representing patients and caregivers were asked to perform simulated injections on both thigh and abdomen.
- Visually and functionally representative models of the devices were used.
- The study included naive and experienced users in equal parts and included some participants with physical impairments.

Result

The PiccoJect™ design was extremely well received. 92% of the participants selected PiccoJect™ when asked for their overall preference.

The wrap-around drug window and the non-round cross-section were particularly appreciated by participants:

- 83% of the participants rated the visibility of the PiccoJect™ drug window as superior.
- 75% of the participants reported that they preferred the wide cross-section of PiccoJect™.

Participants rated the 1 ml variant of PiccoJect™ 28% higher than the comparable marketed device when asked to indicate their willingness to adopt (based on a 7-point Likert scale).

92 %
overall preference
for PiccoJect™

83 %
rated drug window
visibility as
superior

75 %
preferred the wide
cross-section

28 %
higher willingness
to adopt



Haselmeier™ GmbH

Vaihinger Straße 48

70567 Stuttgart, Germany

Phone: +49 711 71978-0

Email: customerservice.drugdelivery@medmix.com

For more information visit us at

www.haselmeier.com

www.medmix.swiss



Disclaimer

The products shown in this product brochure are under development and some of them may not yet have been approved for sale under applicable medical device regulations. The content provides general information about these products, their field of application and intended use, and their function, performance and specification are subject to customer specific development and may deviate from those shown herein. All information contained herein is directed at medical device manufacturers and pharmaceutical companies. This information shall not constitute a promotion for use in a way which conflicts with any applicable medical device regulations, nor is it directed at patients or intended to replace professional medical advice.