




MeCH0ta[®] CHO-K1 Cell Line Development Platform

MeCHOta® Introduction

MeCHOta® CHO-K1 cell line is derived from ECACC CHO-K1, adapted to suspension culture in serum-free medium and selected by monoclonal screening. It has global authorization and proven supporting service kit(expression vector, medium), ensuring stable production capacity from laboratory to commercial scale.

MeCHOta® CHO-K1 Pedigree

MeCHOta® CHO-K1 is derived from CHO-K1 (ECACC 85051005), with clear background and clear lineage, meeting the compliance requirements of excellent commercial engineering cell lines.



Superity of MeCHOta®

Stability, compliance (traceability) and scalability (excellent metabolism) are the most important indicators for evaluating commercial engineering cell lines. MeCHOta® CHO-K1 has a clear background, complete traceability and compliance, has passed comprehensive cell testing, has reliable stability, excellent metabolic level and environmental tolerance, meets all standards of engineering cell lines, and is an excellent commercial project cell line.





MeCH0ta[®] Host Cell Background Traceability



MeCH0ta[®] CHO-K1 has a complete traceability of the whole process from original purchase to cGMP banking and testing, ensuring the compliance of cell line for IND and BLA.

Excellent Growth Characteristics of MeCH0ta[®] Host Cells

MeCH0ta[®] CHO-K1 has uniform cell shape and size (~15 μm), stable growth, short doubling time (~18h), and excellent growth characteristics. It is an ideal cell line for commercial large-scale production.



Complete Detection Report of MeCH0ta[®] Host Cells


The MCB and WCB established by MeCH0ta[®] CHO-K1 under cGMP conditions have been tested for cell identification, sterility, mycoplasma, endogenous and adventitious viruses, and retroviruses, ensuring production compliance and safety under the GMP system.

Sponsor	Thousand Oaks Biologics Inc.			
Sample name	MeCHOTa® CHO-K1 Fus-MCB			
Lot No.	*****			
Start and end date	2023-10-17-2023-12-15			
Test rationale	ChP/USP/CHQ5A(R1)/CHQ5D/US FDA PTC (1997)/US FDA PTC (1993)/9CFR			
Test item	Test method	Test code	Reference standard	Test result ¹
Cell identification	Cytochrome oxidase I (CO1) detection	006002	Source from China Hamster cells	Source from China Hamster cells
	Membrane filtration (Method suitability study)	007004	Report	Method suitable
Test for bacterial and fungi	Membrane filtration	007004	Negative	Negative
	Pouring method (Method suitability study)	008001	Report	Method suitable
Mycobacterium test (Declare for China)	Pouring method	008001	Negative	Negative

Test item	Test method	Test code	Reference standard	Test result ¹
Detection of murine minute virus	Cytopathic observation method (NB324K)*	016002	Negative	Negative
	Quantitative real-time PCR*	016003	Negative	Negative
	Cytopathic effect (BT, MDBK, Vro, MRC-5)	017002	Negative	Negative
	Hemadsorption (BT, MDBK, Vero, MRC-5)	017002	Negative	Negative
Detection of bovine virus	Fluorescent antibody detection:			
	Bovine viral diarrhea virus			
	Bovine parvovirus			
	Bovine adenovirus type 3			
	Reovirus type 3			
	Bovine parainfluenza virus type 3	017002	Negative	Negative
	Bluetongue virus			
Bovine respiratory syncytia virus				
Rabies virus				


MeCHOTa® Perfect matching of host cells

MeCHOTa® CHO-K1 has a proven supporting services kit, including the expression vector system (pTOB series) with independent intellectual property rights, and the corresponding high-expression medium, to help customers quickly advance projects from DNA sequences to shelf products.





Reliable Stability of MeCHOTa® Host Cell

MeCHOTa® CHO-K1 exhibits excellent stability in multiple dimensions such as copy number of exogenous genes, cDNA sequence, passaging, expression level, and expression quality, and no sequence variants are produced, meeting the requirements of commercial continuous production.




Cell Growth of Cell Stability Study



Titer of Cell Stability Study

High expression of MeCHOTa[®] CHO-K1 host cells





Cell lines and types of culture medium

Molecule Code (based on amino acid sequence)	Molecule Type	Cell Line
A	mAb-IgG1	MeCHOTa [®] CHO-K1
B	mAb-IgG1	MeCHOTa [®] CHO-K1
C	mAb-IgG1	MeCHOTa [®] CHO-K1
D	mAb-IgG1	MeCHOTa [®] CHO-K1
E	BsAb-IgG1	MeCHOTa [®] CHO-K1
F	mAb-IgG1	MeCHOTa [®] CHO-K1
G	mAb-IgG4	MeCHOTa [®] CHO-K1
H	mAb-IgG4	MeCHOTa [®] CHO-K1
I	IgE-λ	MeCHOTa [®] CHO-K1
J	mAb-IgG1	MeCHOTa [®] CHO-K1
K	mAb-IgG1	MeCHOTa [®] CHO-K1

Note: MeCHOTa[®] CHO-K1 is derived from ECACC CHO-K1.

Titer and analysis of molecules



Our MeCHOTa[®] CHO-K1 host cells demonstrate high-yield production, achieving high titers ranging across 11 molecules - including 9 monoclonal antibodies, 1 bispecific antibody (8.1 g/L), and 1 IgE-λ molecule.


Notably, molecules D (9.2 g/L) and K (8.9 g/L) showcase exceptional yield performance.

Molecule	Molecule type	Qp (pg/cell/d)	SEC			nrCE-SDS			rCE-SDS
			HMW (%)	Main (%)	LMW (%)	LMW (%)	Main (%)	HMW (%)	LC+HC (%)
C	mAb-IgG1	52	1.8	98.1	0.0	5.7	94.3	/	98.9
D	mAb-IgG1	56	0.0	99.3	0.7	4.0	96.0	/	99.0
G	mAb-IgG4	44	3.7	96.1	0.2	4.3	95.7	/	98.1
H	mAb-IgG4	68	3.9	93.5	2.5	3.9	96.1	/	/
E	mAb-IgG1	32	1.3	98.5	0.2	6.1	93.9	/	98.1

Molecule	Molecule type	N-Glycan (%)										
		G0-GN	G0F-GN	G0	G0F	Man5	G1Fa	G1Fb	G2F	G2S1F	G2S2F	Others
C	mAb-IgG1	1.5	3.5	6.4	54.5	4.8	16.5	6.1	2.0	/	/	4.7
D	mAb-IgG1	1.9	2.9	6.9	40.8	6.0	23.7	6.8	4.1	/	/	6.9
G	mAb-IgG4	0.5	1.1	3.9	71.6	1.3	8.0	7.9	2.1	/	/	3.6
H	mAb-IgG4	1.7	3.5	6.4	48.3	5.1	11.5	11.2	3.9	/	/	8.4
E	mAb-IgG1	0.3	0.6	3.5	51.9	0.8	19.9	8.7	6.5	1.8	1.5	4.5

CHO-K1 Sub-Licensing Process in 4 Steps

Transparent, Efficient, and Risk-Free



(Note: ECACC = European Collection of Authenticated Cell Cultures, a recognized world-leading cell bank)

License from ECACC

Three Key Guarantees

License Reference Number: BE20/123

Licensee's Confirmation:

The Licensee confirms that it has read and understood UK Health Security Agency's Terms & Conditions for Supply of Goods and/or Services (as annexed/printed overleaf).

Authorised signature:	Date: 9th October 2024
On behalf of UKHSA	
Name in block capitals: Judith Parker	In the capacity of: Head of Business Development
Authorised signature:	Date: 2024.09.13
On behalf of the Licensee	
Name in block capitals: SHIN	In the capacity of: president

- Global Authorization:** Worldwide commercialization rights.
- One-Time Fee:** No annual fees, royalties, or hidden costs.
- Termination Protection:** Existing sub-licenses are unaffected by future agreement changes.

Thousand Oaks Biologics Inc.



Building B1, 100 Dongting Lake Road, Linjiang Town, Haimen District, Nantong, Jiangsu Province, China
Building 5, 1215 Zhengjia Road, Lin-gang Special Area, China (Shanghai) Pilot Free Trade Zone
Room 206-207.2nd Floor, No.1043 Harley Road, China(Shanghai) Pilot Free Trade Zone

Email: services@tobiopharm.com Web: www.tobiopharm.com/en