

THE 100% NO TOUCH
TITANIUM HOLDER

THE HOLDER FOR A KLEAN DENTAL SURGERY



Tilt, Hold & Pick safely



PROTECTIVE: The implant is held thanks to a bespoke titanium ring and its own covering/healing screw.



NO TOUCH feature: Surgeons / Scrub nurses can remove directly the implant from the holder without touching it, thus avoiding any contamination.



SMART: Compatible with the «Clean Implant» label. Suitable for right/left handed use.



CUSTOMIZABLE: Can be packaged in pouches, blisters and SoBliss®.



ADJUSTABLE: SoKlean® is able to receive several sizes of implants thanks to the unidirectional slider and our wide range of colored rings.

SoKlean® sterile packaging

Pouch



Blister



so bliss[®]
FLAT



User guide



TILT
LEFT OR RIGHT



HOLD
THE BODY




PICK
SAFELY



PUSH
THE SLIDER



GET
THE SCREW



Implant size

Dmax: 5,5 Lmax: 18
Dmin: 3 Lmin: 6

Ring color



Other colors available upon request



www.selenium-medical.com | contact@selenium-medical.com

Selenium Medical: 40 rue Chef de Baie, 17000 La Rochelle FRANCE | Tel: +33 (0)5 46 44 40 82

Selenium America: 6201 Fairview Road, Suite 200, 28210, Charlotte NC USA | Tel: +1 980 215 8506

