

BD HypakTM for Biotech

Glass Pre-fillable Syringe

Cover your bases in biotech

Grounded on decades of BD Hypak™ experience, BD Hypak™ for Biotech has been developed to address the key compatibility challenges for biotech drugs and chronic treatments.

A reliable Pre-fillable Syringe from a trusted partner

More than 75% of major injectable biotech drugs are delivered in BD Pre-fillable Syringe solution¹

The BD Hypak™ for Biotech is a reliable and proven Sterile, Clean and ready to Fill (BD SCF™) syringe barrel container featuring key compatibility factors for biotech drug delivery:



| BD SCF™* | 1 mL Long |
|-----------------------|-------------------------------|
| Flange | Cut or Small Round |
| Needle Cap type | Rigid Needle Shield |
| Needle diameter | 27G – 29G |
| Needle Wall thickness | Regular, Thin or Special thin |
| Needle length | 1/2inch |
| Needle type | 5 bevels |

BD Medical

Pharmaceutical Systems **United States** 1 Becton Drive Franklin Lakes, NJ 07417 +1 800 225 3310

Europe

11 rue Aristide-Bergès 38800 Le Pont-de-Claix

Phone: +33 4 76 68 36 36 - Fax: +33 4 76 68 35 05 Becton Dickinson France S.A.S - Share capital: 64 719 915 €

RCS Grenoble B 056 501 711

and gliding specifications foundational for compatibility with secondary devices 3,4



Disposable Autoinjector

Needle guard

References and notes

- * SCF: Sterile, Clean and ready to Fill (BD SCF^{TM}).
- 1. Internal analysis based on revenues projections in the scope of 15 therapeutic classes in the chronic conditions, BD internal sales (Vendavo 2017) and BD internal market model (InThought 2017).
- 2. Validation Summary Low Tungsten BD Hypak™ for Biotech Prefillable Syringe Barrels [Internal Study] Pont-de-Claix. FR: Becton, Dickinson and Company; 2012.
- BD Hypak™ for Biotech customer specifications.
- 4. BD Hypak™ for Biotech 1mL Long Prefillage Syringe [Internal Study] Pont-de-Claix, FR: Becton, Dickinson and Company; 2013.



bd.com