



	SAN-Light	SAN-Light DV	SAN-Light DV Pro	SAN-L	SAN-L DV	SAN-L DV Pro	SAN-P	SAN-In
<b>Features*</b>	Passive safety needle	Passive safety needle with vial adaptor	Passive safety needle with integrated vial transfer system	Passive automatic safety needle	Passive Automatic safety needle with vial adaptor	Passive Automatic safety needle with integrated vial transfer system	Automatic safety syringe for staked needle prefilled syringe	Prefilled syringe with integrated safety needle
Primary Container Compatibility	Luer-Lock Syringe	Vial	Vial	Luer-Lock Syringe	Vial	Vial	Staked needle prefilled syringe	COP / COC prefilled syringe
Drug Formulation	Liquid / Lyophilized***	Liquid / Lyophilized	Liquid / Lyophilized	Liquid / Lyophilized	Liquid / Lyophilized	Liquid / Lyophilized	Liquid	Liquid
Possible Customization	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
End User	HCP / Caregiver	HCP / Caregiver	HCP / Caregiver	Self-Injection / Caregiver	Self-Injection / Caregiver	Self-Injection / Caregiver	Self-Injection / Caregiver	HCP / Caregiver
Route of Administration - SC/IM/ID**	SC / IM / ID	SC / IM / ID	SC / IM / ID	SC / IM / ID	SC / IM / ID	SC / IM / ID	SC / IM / ID	SC / IM / ID
Automatic Passive Needle Shielding After Injection	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Automatic Needle Insertion	No	No	No	Yes	Yes	Yes	Yes	No
Manual control of injection speed	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Hidden Needle	Partially	Partially	Partially	Yes	Yes	Yes	Yes	Partially

\* Each device can be customized according to customer's requirements; e.g., needle gauge and length, industrial design, primary container, etc.  
 \*\* SC - subcutaneous; IM - intramuscular; ID - intradermal  
 \*\*\* Usage with Dual-Chamber Syringes