
Centurion Self-Injection Systems:

2-Step Auto-Injector

- ✓ Developed and tested.
- ✓ Looking collaboration for scale-up and CE processes



Flexible delivery volumes:0.1 - 1.0 mL standard 1 mL PFS

2 step easy operation:Remove orange capPush against injection site

Automatic needle insertion and injection

Passive or active needle guard options



Centurion Self-Injection Systems:

Centruion PFS Support System



- Flexible delivery volumes:0.1 1.0 mL long 1 mL PFS
- 2 Step easy operation:
 - Remove blue cap
 - Push against injection site
 - Inject the solution

Needle hidden system before and after injection

Small size and reduced cost

Contact: Associate Professor Adem Sahin ademsahin@centurion.com.tr







