



ONE2SAMPLE™

SIMPLE, SECURE AND CONTROLLED SAMPLING



THE SANISURE ADVANTAGE

Your Go-To Partner for Flow Path and Closed-System Solutions

Created from the merger of six companies leading the market in their fields of expertise, SaniSure is unmatched in our ability to provide you with a uniquely diversified and vertically integrated portfolio of highly innovative proprietary solutions, custom assemblies, and engineered components.

We apply our decades-deep product and application expertise in a consultative manner, working with you to resolve your most challenging flow path needs, supported by our industry-leading, in-house design and manufacturing capabilities.

While larger companies cannot offer the flexibility, and smaller companies lack the depth of capabilities, we stand apart as the best of both worlds—unparalleled flexibility combined with end-to-end capabilities and deep vertical integration to help you de-risk, grow and innovate.



ONE2SAMPLE™

Draw up to 10 samples with a single connection using configurable assemblies with a variety of container types and combinations. One2Sample is the only solution that streamlines workflows with secure, flexible and efficient simultaneous or sequential sampling.



SECURE. FAST. EASILY CONFIGURABLE.

Secure

Sterile, Representative Samples

- Single-connection filled samples ensure each sample is identical and representative.

Reduced Process Risk

- A single, simple and reliable sterile connection reduces the risk of process contamination by up to 10X versus single sampling.
- One pre-sterilized assembly of up to 10 samples reduces sterility assurance management.
- One2Sample hub and Cap2V8® bottles and caps to sampling vessel offers a completely closed system.
- A single Aseptiquik connection enables filling 10 pre-sterilized containers without the use of a biosafety cabinet.
- Pinch pipes offer fast and reliable aseptic disconnection.

Fast

10X Time Savings

- One2Sample fills 10 samples with a single connection, reducing time spent by up to 90%.



Easily Configurable

A wide variety and mix of sampling containers means every One2Sample custom assembly can be configured based on the specific unit operations' sampling needs.

1

Every sample is in the right container for the analytical method.

2

Each sample can be filled directly with no additional handling or sample prep.

3

Simplified handling reduces opportunities for sampling errors.

SAMPLING SOLUTIONS FOR LIFE-SAVING BIOPHARMACEUTICALS

As a trusted partner for biopharmaceutical companies worldwide, SaniSure offers secure, sterile, easily configurable sampling with countless sampling container options to perfectly suit your bioprocessing needs.

Biopharma Processes

Upstream	Media Preparation: Sampling for Physicochemical and Sterility
Harvest	Harvest Pool: Sampling for Physicochemical, Impurities, and Bioburden
Downstream	Chromatography Pooled Fractions: Sampling for Physicochemical, Identify, Potency, Impurities, and Sterility Filtration (Protein A and Viral Inactivation): Sampling for Physicochemical, Impurities, Identify, and Sterility UDFD (Start and Finish): Sampling for Physicochemical, Impurities, Identify, and Sterility Intermediates: Sampling to Ensure Sterility and Stability
Final Fill	Final Formulation: Sampling for Physicochemical, Potency, Impurities, Identify, Sterility, and Stability Samples Only one vent filter to test for integrity if required by your QA.
Custom Configurations	Assemblies are configured to your specific, single-use process needs. All components are readily available with compliance data.

PROCESS CONNECTION OPTIONS

One2Sample assemblies are easily configurable to include a vast variety of connectors that meet your specific sampling needs. They also offer fast and reliable aseptic disconnection with pinch pipes and thermal seals on CellGyn® tubing.



Weldable Tube

Process Connection Tubing Sizes Available:

- 1/8" x 1/4"
- 1/4" x 3/8"
- 3/8" x 5/8"



Quick Connect

Process Connection Tubing Sizes Available:

- 1/8", 1/4", 3/8" ID



Sterile Connectors

Process Connection Tubing Sizes Available:

- 1/8", 1/4", 3/8" ID
- 3/8", 1/2" ID






Sanitary Flange

Process Connection Tubing Sizes Available:

- 3/4" TC: 1/8" – 1/2" ID
- 1" – 1/2" TC: 1/4" – 1" ID

SINGLE VENT HUB

Easily adaptable to your vessel size with multiple fluid line options, we significantly reduce the effort and time required for vent filter integrity testing, minimizing this effort by a factor of 10 with our One Vent Hub solution.



SAMPLE CONTAINERS

A wide variety and mix of sampling containers means every sample is in the right container for the analytical method.

Conical
Tubes
Starting at 5 mL



Media
Bottles
Starting at 60 mL



Port
Bags
Starting at 50 mL



	SWABABLE LUER	PP CONICAL TUBE	PC & PET PHARMATAINER™	PETG MEDIA BOTTLE	2 PORT BAG
SIZES AVAILABLE	Female Luer	5 mL, 15 mL, 50 mL	50 mL, 125 mL, 250 mL	60 mL	50 mL, 150 mL, 250 mL
VENT	N/A	50 mm Shared Vent Filter			N/A
ASEPTIC DISCONNECTION	Thermal Seal on CellGyn or Pinch Pipe				
SPECIFIC APPLICATIONS	Syringe access point for pH, DO, DCO2	General Use	PET: Stability, Bioburden, Retain PC: Bioburden, Endotoxin	Customer Specific Request	General Use
KEY FEATURE	Flexible, Accurate Volume	Low Cost, Accurate Volume	Clean, Clear, Low Endotoxin Binding	Clear	Oversized

INTEGRATE SEAMLESSLY

Proprietary Solutions



Mixed4Sure™



One2Fill™



Fill4Sure™



RC-40™



Custom Assemblies

Engineered Components



Bottles and Stability Vials



Cap2V8® Bottle and Cap Assemblies



Diversified Tubing Portfolio



aSURE™ Instrument Tees



Closed System Vessels



2D and 3D Bags



Validated Library of Fittings, Clamps and Other Components



Precision Filling Needles



sanisure.com

Americas

1020 Flynn Road
Camarillo, CA 93012
United States

+1 800 726 4835

SalesSupport-US@sanisure.com

Asia-Pacific

Suite 708, German Center
88 Keyuan Road
Pudong District, Shanghai, PRC

+86 21 68771098

SalesSupport-CN@sanisure.com

Europe, Middle East, & Africa

ZAE Robert Steichen
16, rue Héierchen
L4940 Bascharage
Grand Duché du Luxembourg

+352 26 33 731

SalesSupport-EU@sanisure.com