


OXYCAPT™ Multilayer Plastic Vial

Multilayer Structure




- Excellent Oxygen Barrier
- High Water Vapor Barrier
- Very Low Extractables
- Low Protein Adsorption
- Excellent Ultraviolet Barrier
- High Break Resistance
- High pH Stability
- Gamma-sterilized Vial
- Customizable
- For Biologics & Regenerative Medicine




OXYCAPT™ Multilayer Plastic Vial


Oxygen Barrier (23℃・75%RH)



Ultraviolet Barrier



Inorganic Extractables (ICP-MS)



Product Line-up

ISO Vial	Length (mm)	Width (mm)	OD of Flange (mm)	ID of Flange (mm)	Option
2R (2mL)	35	16	13	7	Bulk or RTU
6R (6mL)	40	22	20	12.6	Bulk or RTU
10R (10mL)	45	24	20	12.6	Bulk or RTU
20R (20mL)	55	30	20	12.6	Bulk or RTU

*These data are general test results by standard test methods and not guaranteed values.



Mitsubishi Gas Chemical Company, Inc.
New Business Development Department
Mitsubishi Building, 5-2 Marunouchi 2, Chiyoda-ku, Tokyo 100-8324, Japan
Phone: +81 3 3283 4913
<https://www.mgc.co.jp/eng/products/abd/oxycapt.html>



(Issued in June 2022)