



About ACT



ACT is a supplier of a variety of premium empty capsules, offering **H-Caps®(HPMC)**, **H-Caps®DR**, **Gelatin**, **Pullulan**, and custom-made capsules to meet the needs of the nutraceutical and pharmaceutical industries. With a commitment to quality and innovation, we continuously strive to deliver the highest standard of capsules available.

Our capsules are available in a range of sizes, including **#00EL**, **#00**, **#0EL**, **#0**, **#1**, **#2**, **#3**, **#4**, and **#5**, ensuring flexibility to meet diverse formulation requirements.

Add: Eckle 7, 77704 Oberkirch, Germany

Tel: [+49 1782331795](tel:+491782331795)

Amtsgericht Freiburg i.Br. HRB 730767

Ust-IdNr.: DE369149022/ BE1016469433

Steuer Nr.: 14002/65220

Sparkasse Offenburg/Ortenau

IBAN:DE98664500500004987882

BIC: SOLADES1OFG

EORI-Nr. DE443822072544895

ACT & GDL



ACT serves as the sole European distributor of GDL Capsules, bringing GDL's premium capsule products to the local market.

As the European extension of GDL Capsules, ACT bridges the gap between GDL's expertise and local customers, ensuring seamless access to top-quality hard capsules.



Our commitment to outstanding service is reinforced by extensive industry knowledge, guaranteeing reliability and excellence in every product we provide.



Overview of Branches



(1982) GDL Capsule (Taiwan) Co., Ltd.

GDL was originally founded and established in Taiwan, where the company officially began manufacturing empty two-piece hard capsules in 1982. Utilizing self-developed machines and technologies, GDL ensured a high standard of quality from the start.

(1992) GDL Capsule (Shanghai) Co., Ltd.

In 1992, Shanghai GDL factory was established to broaden its reach across mainland China and to cater the fast-growing needs of hard capsules in the pharmaceutical industry.

(2003) GDL Capsule Phils., Inc. (Philippines)

GDL extended its sales network to Southeast Asia by establishing an office in Manila in 2003, enhancing its ability to support regional customers.

(2008) Biogel Capsules Corporation (North America)

Continuing its international expansion, GDL set up Biogel Capsules Corporation in 2008 as a dedicated sales office in the United States, providing local support to North American clients.

(2019) GDL Pharmaceuticals (Guizhou) Co., Ltd.

This new subsidiary represents the latest technology of GDL with streamline production process and Hard Capsule Machines from sister company, GDW. Operations began in year 2021 and an expansion of capacity is currently carrying out to satisfy the growing demand of vegetarian HPMC and Pullulan Capsules.

(2024/2025) ACT Capsules GmbH (Germany)

ACT Capsules GmbH founded in Germany in March 2024 is to serve the European customers more effectively. Official sales are slated to begin in January 2025.

GMP complied and ISO certified



All Capsules are Free from ETO and Radiation Treatments, with Production Processes Strictly Controlled and ISO 9001:2015 Certified

US DMF No:
Gelatin capsule: 12611 HPMC capsule: 034924

Certification



CERTIFICATE

TÜV SÜD Management Service GmbH (043)
certifies that, having conducted an audit
the company

GUIZHOU GUANG DE LI PHARMACEUTICALS CO., LTD.
Daping Industrial Park, Economic and Technological Development Zone
Qiannan Buyi and Miao Autonomous Prefecture,
Guizhou Province, P.R. China
Post Code: 558000
Site code: 10005830

meets the requirements set out in the

**GLOBAL STANDARD FOOD SAFETY
Issue 9: August 2022.**

Scope: Production (moulding, drying) of
hydroxypropyl methylcellulose (HPMC) capsules, gelatin vacant capsules,
pullulan vacant capsules and plants vacant capsules packed in PE plastic bag
Product Categories: 15
Achieved Grade: A
Audit program: Announced

Date of Audit: 2025-07-16 - 2025-07-18
Re-audit due date: 2026-06-27 - 2026-07-25

Certificate Registration No.: 12 020 64509 TMS
Plant No.: 14300013
Auditor No.: 211020
Issue date: 2025-08-08
Expiry date: 2026-09-05

If you would like to feedback comments on the BRCS Standard
of the audit process directly to BRCS, please contact audit@brcs.com or via
the BRCS Reporting Platform at <https://reporting.platform.brscs.com>
To verify certificate validity, please visit <https://www.brscs.com/verify>



PRODUCTS



Gelatin Capsules

- ✓ No BSE/TSE risks
- ✓ No preservatives
- ✓ Halal Certified
- ✓ Excellent filling yield rate
- ✓ Tighter weight specifications can be provided
- ✓ Easy to swallow

TiO₂-Free capsules are available



H-Caps®



(HPMC Capsules)

- ✓ No BSE/TSE risks
- ✓ Free of preservatives, allergens and starches
- ✓ Non-GMO
- ✓ Avoiding risks of cross-linking reactions of gelatin
- ✓ Suitable for water sensitive formulation
- ✓ Non-Animal derived
- ✓ Halal Certified
- ✓ Kosher Certified
- ✓ Reliable quality and filling performance

TiO₂-Free capsules are available

Presented by GDL®

HPMC Capsules



Prulan Caps® (Pullulan capsules)

- ✓ Pure plant source
- ✓ Non-GMO
- ✓ None BSE/TSE risks
- ✓ Low oxygen permeability
- ✓ Free of preservatives, gluten and starch
- ✓ Suitable for organic products
- ✓ Halal Certified
- ✓ Kosher Certified
- ✓ Excellent choice for green consumers
- ✓ Organic certification available upon request.

TiO₂-Free capsules are available



Prulan Cap® (Pullulan capsules)

Excellent choice for Green Consumers

H-Caps® DR Prulan Caps® DR (Delayed-Release Capsules)

- ✓ Delayed release
- ✓ Resistant to gastric acids
- ✓ Resistant to cross-linking
- ✓ Low moisture content
- ✓ Stability over a wide range of temperatures and humidity
- ✓ Optimize the performance of probiotic supplements

TiO₂-Free capsules are available

Capsules for Probiotics & Enzymes



Special Capsules

- SLS Free Capsules
- Pearl Color Capsules
- Fish Gelatin Capsules
- Gelatin Capsules FD (Fast Disintegrating)
- H-Cap®LWC (Low Water Content Capsules)
- H-Cap®DPI (Dry Powder Inhalation Capsules)




Liquid Filled Hard Capsules



Liquid-Filled Capsule



Capsule in-capsule



Beadlets + Liquid

LIQUID CAPSULE FILLING SOLUTIONS

CUSTOMIZED OEM SERVICE FOR DIETARY SUPPLEMENTS

Featured Formulas

1. Deep Sea Fish Oil (DHA + EPA)
(High-Potency)

2. Fish Oil + CoQ10
(Cardiovascular Support Formula)

3. MCT Oil + Micro-Encapsulated Vitamin C & Beta-Carotene
(Advanced Energy & Antioxidant Blend)



Customized Service

Printing

There are options of axial, radial, spin and two-color printing for branding identification.



Rectified Axial Print



Axial Print



Rotational Print


Color

We can customize a wide range of capsule colors for your needs.

You can also choose from our existing capsule colors.



Capsule Sizes



Size	Length(mm)				Weight(mg)	
	Cap	Body	Pre-Locked	Locked	Gelatin Capsules	HPMC Capsules
#00EL	12.20±0.40	21.60±0.40	26.95±0.85	24.65±0.40	134±9	137±10
#00	11.70±0.40	20.10±0.40	25.80±0.85	23.40±0.40	125±8	124±9
#0EL	11.70±0.40	20.20±0.40	25.00±0.85	22.80±0.40	107±8	105±9
#0	11.00±0.40	18.40±0.40	23.75±0.85	21.30±0.40	100±7	96±8
#1	9.80±0.40	16.50±0.40	21.35±0.85	19.20±0.40	78±6	76±7
#2	9.00±0.40	15.30±0.40	19.90±0.65	17.60±0.40	63±5	61±6
#3	8.20±0.40	13.60±0.40	17.75±0.75	15.90±0.40	50±4	49±5
#4	7.30±0.40	12.35±0.40	16.50±0.50	14.50±0.40	39±3	38±4
#5	6.20±0.40	9.30±0.40	12.90±0.30	11.00±0.40	28±3	

Packaging Method

- Pharmaceutical grade PE bag wrapped with aluminum foil bag contained in 5 layer corrugated carton

Storage Conditions

- Warehouse temperature: 15-25°C, Relative humidity: 35-65%
- Cartons should be kept away from direct heat source or sun light
- Keep PE bags sealed at all time with capsules inside
- It is not recommended to use capsules at relative humidity which is higher than 55%

Loading Method

- Carton Dimension: 58.5x39.5x73 (cm)
- 15 Cartons/Pallet

Size	q'ty(pcs)
#00EL	70,000
#00	75,000
#0EL	100,000
#0	100,000
#1	130,000
#2	180,000
#3	230,000
#4	300,000
#5	400,000

Container Type	Loading Method	Capsules/Per Carton
20 feet reefer container	Pallet	108
	No Pallet	147
40 feet reefer container	Pallet	300
	No Pallet	353

Why ACT?

- Assurance of Supply
- High and Consistant Quality
- Excellent Filling Performance
- Optimal Price-Performance Ratio
- Regulatory Compliance
- Formulation Assistance
- Technical Service





Contact us

ACT Capsules GmbH

Dußlinger Weg 10, 72072 Tübingen, Germany

Mobile: +49 1782331795

Email: info@actcap.de

Website: www.actcap.de



Thank you
Vielen Dank
Muchas Gracias