



# Maggie<sup>®</sup>

Two Step Cartridge-Based Solution

5.0  
mL



Volume 

Viscosity 

Device Type	Disposable autoinjector
Primary Container	5.0 mL cartridge with pre-installed needle (NIT <sup>®</sup> )
Dosage Type	Fixed dose   Optional automated priming and titration
Route of Administration	Subcutaneous
Activation Method	Needle cover
Injection Feedback	Continuous audible clicks during injection, viewing window

## Handling Instructions



1  
Unscrew cap until needle cover automatically extends\*.  
Then remove cap.



2  
Press and hold against injection site to activate.



+  
Discard in sharps container.

\* For suspension drugs that require a priming step – hold device with cap upright during cap removal.

## Key Features

- ✓ Simple 2-step operation
- ✓ Needle cover activation
- ✓ Continuous audible and visual feedback
- ✓ One-handed operation
- ✓ Permanently hidden needle
- ✓ Compact design
- ✓ Customizable needle length and gauge
- ✓ Passive locking needle cover
- ✓ Customized to client specifications



The SHL Medical logo and Maggie 5.0 are registered and/or non-registered trademarks of SHL Medical AG or its affiliates, in the United States, Europe, Taiwan, and other jurisdictions. Final product specifications are subject to change based on customer requirements and SHL's offerings. Contact [info@shl-medical.com](mailto:info@shl-medical.com) for more information.

[www.shl-medical.com](http://www.shl-medical.com)