



Molly[®]
Modular Platform
Technology



Volume 

Viscosity 

Device Type	Disposable autoinjector
Primary Container	2.25 mL pre-filled syringe with staked needle
Dosage Type	Fixed Full contents
Route of Administration	Subcutaneous
Activation Method	Needle cover
Injection Feedback	Audible click at beginning and toward end of injection, viewing window

Handling Instructions



Remove cap.



Press and hold against injection site to activate.



Discard in sharps container.

Key Features

- ✓ Simple 2-step operation
- ✓ Needle cover activation
- ✓ Audible, visual, and tactile feedback
- ✓ One-handed operation
- ✓ Larger volume
- ✓ Permanently hidden needle
- ✓ Ergonomic cap removal
- ✓ Robust anti-roll design



The SHL Medical logo and Molly 2.25 are registered and/or non-registered trademarks of SHL Medical AG or its affiliates, in the United States, Europe, Taiwan, and other jurisdictions. Customized designs are for reference only. Final product specifications are subject to change based on customer requirements and SHL's offerings. Contact info@shl-medical.com for more information.

www.shl-medical.com