



**Molly**®  
Modular Platform  
Technology



Volume 

Viscosity 

<b>Device Type</b>	Disposable autoinjector
<b>Primary Container</b>	1.0 mL long pre-filled syringe with staked needle
<b>Dosage Type</b>	Fixed   Full contents
<b>Route of Administration</b>	Subcutaneous
<b>Activation Method</b>	Needle cover
<b>Injection Feedback</b>	Audible click at beginning and toward end of injection, viewing window

## Handling Instructions



1  
Remove cap.



2  
Press and hold against injection site to activate.



+

Discard in sharps container.

## Key Features

- ✓ Simple 2-step operation
- ✓ Needle cover activation
- ✓ Audible, visual, and tactile feedback
- ✓ One-handed operation
- ✓ Ultra-compact design
- ✓ Permanently hidden needle
- ✓ Passive locking needle cover
- ✓ Preconfigured or customized to client specifications