



BIOCATALYTIC CASCADES

TURN THIS...



...INTO THIS

FOR MAKING YOUR APIs!



ONE STEP

~ 20°C (room temperature)
Atmospheric pressure

Find out
more here!



OUR EXPERIENCE = YOUR POSSIBILITIES

P450s/ Peroxygenases	α-Oxoglutarate dep dioxygenases (2OGDs)	Alcohol oxidases	Amino acid deaminases	Amino acid dehydrogenases	Reductive aminases (RedAm), Imine reductases (IRED)	Laccases
Carboxylic acid reductases (CARs)	Monoamine oxidases (MAOs)	Transaminases (TA)	Tyrosine-phenol lyase (TPL)	Tyrosine- ammonia lyase	Hydratases	Aldolases
Carboligases	Hydroxynitrile lyase (HNL)	ATP-dependent ligases X = N, O, S	Phenol (De)carboxylases	Fatty acid photo- decarboxylases	Ketoacid decarboxylases	Thiolases
Alcohol dehydrogenases	Halohydrin dehalogenases X = N, O, S	Epoxidases	Prenyl- transferases	Ene-reductases	Sulfatases	Kinases / Phosphatases