



BIOCATALYTIC CASCADES

TURN THIS...



...INTO THIS

FOR MAKING YOUR APIs!



ONE STEP

~ 20°C (room temperature)
Atmospheric pressure

Find out
more here!



OUR EXPERIENCE = YOUR POSSIBILITIES

<p>P450s/ Peroxygenases</p>	<p>α-Oxoglutarate dep dioxygenases (2OGDs)</p>	<p>Alcohol oxidases</p>	<p>Amino acid deaminases</p>	<p>Amino acid dehydrogenases</p>	<p>Reductive aminases (RedAm), Imine reductases (IRED)</p>	<p>Laccases</p>
<p>Carboxylic acid reductases (CARs)</p>	<p>Monoamine oxidases (MAOs)</p>	<p>Transaminases (TA)</p>	<p>Tyrosine-phenol lyase (TPL)</p>	<p>Tyrosine- ammonia lyase</p>	<p>Hydratases</p>	<p>Aldolases</p>
<p>Carbolligases</p>	<p>Hydroxynitrile lyase (HNL)</p>	<p>ATP-dependent ligases X = N, O, S</p>	<p>Phenol (De)carboxylases</p>	<p>Fatty acid photo- decarboxylases</p>	<p>Ketoacid decarboxylases</p>	<p>Thiolases</p>
<p>Alcohol dehydrogenases</p>	<p>Halohydrin dehalogenases X = N, O, S</p>	<p>Epoxidases</p>	<p>Prenyl- transferases</p>	<p>Ene-reductases</p>	<p>Sulfatases</p>	<p>Kinases / Phosphatases</p>