

One2Fill™ 10x Simultaneous Fill



One2Fill™ (Simultaneous Fill) is a range of standard manifolds with bottles for aliquoting or sampling in a variety of upstream and downstream applications, including media or buffer preparation, and other filling applications.

- One vent hub = one vent filter test per assembly
- One feed line to simplify your filling process and reduce the risk of error
- Easily handleable manifold compact trays available to secure your filling operations
- Filling accuracy of $\leq 3\%$

Features

Benefits

One Vent Hub

- Innovative solution patented by SaniSure[®].
- Only one vent filter per manifold ensures a reduction in filter integrity testing and lowers overall costs.
- Robust connection with hub system: cable tie free.

PharmaTainer™

- Available in PC or PET depending on application needs.
- 125ml, 250ml, 500ml, 1L and 2L.

Cap2V8[®] + Cellgyn[®]

- One piece molded cap system (cable tie free) reduces risk of contamination and leakages.
- E&L and compatibility studies available.

Pinch Pipe

- Aseptic disconnection.
- Risk of cross contaminations reduced.

Aseptiquik G from Colder

- Enable quick and easy sterile connections.
- Reliable performance.
- Quick and easy: only one sterile connection needed

Tray

- Robust
- 10 bottles
- Small footprint
- Re-usable

One2Fill® 10x Simultaneous Fill

Packaging

- Double-bagged and labelled in an ISO Class 7 cleanroom

Sterilization


- Gamma radiation compatible (@50kGys)
- Sterility assurance level 10⁻⁶

Materials


- ① 10 bottles 125, 500, 1000 or 2000ml in PET or PC
- ② Cap2V8® 48 Y 2 24" W/ 1/8" ID x 1/4" OD tubing, Cellgyn®
- ③ 3/8" ID X 5/8" OD tubing, Cellgyn® (CGN 375X625)
- ④ 50mm 0.2µm hydrophobic vent filter, gamma-stable (A50V / 17805 AI)
- ⑤ 3/8" HB x Aseptiquik G connector, Polycarbonate (AQG17006)
- ⑥ 10 Port Hub x 0.125" HB Double Kit, Polypro
- ⑦ 1/4" Pinch pipe (NS6 15 05)
- ⑧ Small Pinch clamp, Polypropylene (340TCSP)
- ⑨ LG. Pinch clamp, Polyester (340TCL)
- ⑩ STD. Cable tie, Nylon (BT2S M39)

Applications

- Media and Buffer Aliquoting
- Sampling
- Drug Substance Filling



Trays – reusable (Polypropylene and stainless steel),
Delivered separately from manifold



Drawings on request
E.g. Tray for 10 way One2Fill™ 125 PET or PC – SS96200S

Ordering

Volume	Material	Manifold Part Number	Tray Part Number
125 ml	PET	SS99030S	SS96200S
	PC	SS89995S	
250 ml	PET	SS99031S	SS96205S
	PC	SS89996S	
500 ml	PET	SS99035S	SS96210S
	PC	SS99000S	
1000 ml	PET	SS99040S	SS96220S
	PC	SS99005S	
2000 ml	PET	SS99045S	SS96230S
	PC	SS99010S	

One2Fill™ 10x can be customized according your needs.
Please contact your sales representative.

Americas

1020 Flynn Road
Camarillo, CA 93012
United States

+1-805-920-7264
SalesSupport-US@SaniSure.com

Europe & Asia

ZAE Robert Steichen 16, rue Héierchen
L4940 Bascharage
Grand Duché du Luxembourg

+352-26-33731
SalesSupport-EU@SaniSure.com



SaniSure.com